



# Thermo

ELECTRON CORPORATION

Potentiometric Titration Application Notes

Applications Log # 316

## Overview

N, N-bis (2-hydroxyethyl alkyl amine) concentration was determined using the first derivative titration technique with an Orion pH electrode, and perchloric acid as the titrant. The Orion 960 Autochemistry System determines the endpoint and calculates the concentration of N, N-bis alkyl amine in the sample.

---

<b>Industry</b>	Plastics, Rubber, Polymers				
<b>Species Measured</b>	Alkyl Amine				
<b>Sample</b>	Polyethylene				
<b>Sample Size</b>	5.0g				
<b>Typical Concentration</b>	1662ppm(w)				
<b>Technique</b>	# 6 First Derivative				
<b>Electrode</b>	Ross Sure-Flow Combination pH 8172BN				
<b>Solutions</b>	Deionized water; perchloric acid; glacial acetic acid; n-butyl alcohol; reagent grade methanol; potassium iodide (pellets)				
<b>Sample Prep</b>	Accurately weigh out 5 g of sample and grind to a powder if possible. Add ground sample to an analysis beaker containing 50 mL of n-butyl alcohol and stir thoroughly. After stirring, the sample may be filtered on a suction filtration assemble and the filtrate collected for analysis or the sample may be analyzed as is. This sample preparation may differ depending on user's specification of sample preparation.				
<b>Statistics</b>					
<b># of Trials</b>	5	<b>Mean</b>	1661.5ppm(w)	<b>%CV</b>	1.13
<b>Analysis Time</b>	3.0minute(s)				
<b>Comments</b>	Rinse the electrodes, stirrer, and dispenser probe between measurements with reagent grade methanol.				

## Method Parameters

<b>Sample Volume/Weight</b>	5.1 g	<b>Timed or Stability Readings</b>	5.0 second(s) timed
<b>Constant Increment</b>	0.295 mL	<b>Number of Endpoints</b>	1
<b>Max Titrant Volume</b>	15.00 mL	<b>Desired Units</b>	ppm - w
<b>Molecular weight</b>	349.84 g	<b>Predose</b>	0.098 mL
<b>Prestir</b>	10.0 second(s)	<b>Additional Parameters</b>	
<b>Reaction Ratio</b>	1.00		