



Thermo

ELECTRON CORPORATION

Potentiometric Titration Application Notes Applications Log # 563

Overview The concentration of anionic surfactant was determined by the first derivative technique utilizing Hyamine 1622 as the titrant. The Orion 960 Autochemistry System determines the endpoint and calculates the sample concentrations.

Industry	Chemical
Species Measured	Anionic Surfactant
Sample	Petroleum/Metal Cutting Drawing
Sample Size	4mL
Typical Concentration	8.00%
Technique	# 6 First Derivative
Electrode	Surfactant electrode 9342BN. D-J Ref electrode 900200
Solutions	Ref electrode fill 810007-900002. Hyamine 654201. Triton 654203. 0.005M HCl rinse
Sample Prep	Pipet 4 mL of sample into 150mL beaker. Dilute w/deionized water and then add 2mL triton. Titrate w/ hyamine.

Statistics

# of Trials	3	Mean	7.18	%CV	0.50
Analysis Time	6.9minute(s)				
Comments	Rinse the electrodes, stirrer, and dispenser probe between measurements with deionized water.				