



Thermo

ELECTRON CORPORATION

Potentiometric Titration Application Notes

Applications Log # 676

Overview

The calcium content in a 20 ppm standard solution determined via an indicator titration using the Orion 960 Autochemistry system and an Orion 9629BN cupric combination electrode. The standard solution electrode. The standard solution also contains 2 ppm cupric as the indicator, EGTA is used as the titrant, and the endpoint is determined using a first derivative technique.

Industry	Chemical
Species Measured	Calcium
Sample	20 ppm standard solution of Calci
Sample Size	5.0 ml
Typical Concentration	~20.0 ppm
Technique	# 6 First Derivative
Electrode	Cupric Combination electrode
Solutions	Internal Filling Solution (Cat. 900063); 0.1 M Calcium standard (Cat. 922006); 0.1 M Cupric standard (Cat. 942906); Sodium Tetraborate decahydrate (Na ₂ B ₄ O ₇ 10H ₂ O); 0.0009874 M EGTA; DI water. Thermo Orion 960 Autochemistry System (Cat. 096000);
Sample Prep	Pipette 5 ml standard solution in a 150 ml plastic beaker. Add 45 ml of tetraborate decahydrate pH solution in the above beaker. Adjust pH between 10.5~11.0 by using ammonium hydroxide.

Statistics

# of Trials	5	Mean	21.8 ppm	%CV	0.32
--------------------	---	-------------	----------	------------	------

Analysis Time 10.7 min

Comments Rinse the electrodes, stirrer, and dispenser probe thoroughly between measurements with deionized water.
Bias in result is attributed to 2 ppm cupric present.

Method Parameters

Sample Volume/Weight	5 ml	Timed or Stability Readings	10.0 sec
Constant Increment	0.1 ml	Number of Endpoints	1
Max Titrant Volume	4.0 ml	Desired Units	ppm (v)
Molecular weight	40.0	Predose	1.0 ml
Prestir	3.0sec	Additional Parameters	Total solution volume 50 ml; Titrant concentration 0.0009874 M of EGTA.
Reaction Ratio	1.0		