



Thermo SCIENTIFIC

Potentiometric Titration Application Notes

Applications Log # 641

Overview

Anionic surfactant in shampoo can be titrated with 0.05M Hyamine using the first derivative titration technique, the Orion 960 Autotitrator PLUS, and the Orion Surfactant electrode.

Market	Consumer Products	Species Measured	Anionic Surfactant
Sample	Shampoo	Sample Size	0.5g
		Typical Concentration	< 10% w/w
Technique #	66 First Derivative	Electrode	Surfactant electrode 9342BN. D-J Ref electrode 900200
Solutions	Inner/Outer fill solution 900002, 810007. Hyamine 654201. SLS Standard 654202.		
Solutions preparation:	None required.		
Titrant standardization	Dilute 15.0 mL of SLS standard with about 85 mls of DI water in a beaker. Titrate.		
Sample Prep	Weigh about 0.5 g and record exact weight. Add 100mL deionized water to sample. Sample is now ready to titrate.		
Statistics			
# of Trials	3	Mean	6.70% w/w
		%CV	0.17
		Analysis Time	8.4 minute(s)
Comments	Rinse the electrodes, stirrer, and dispenser probe between measurements with deionized water.		

Method Parameters

Sample Volume/Weight	0.86 g	Timed or Stability Readings	10.0 mV/min stability
Constant Increment	0.250 mL	Number of Endpoints	1
Max Titrant Volume	5.00 mL	Desired Units	% w/w
Molecular weight	370.0	Predose	none
Prestir	1.0 second(s)	Additional Parameters	
Reaction Ratio	1.00		