



# Thermo SCIENTIFIC

*Karl Fischer Titrator*

*Log #* KF1579

**Overview** The moisture in starch powder in gelatin capsules was determined by Karl Fischer titration on the Orion Turbo2 Blending Volumetric Titrator equipped with a high speed blender. Methanol-formamide mixture was used as the solvent.

**Sample** Starch in gelatin capsules

**Industry** PHARMACEUTICALS

**Typical moisture** 8 - 9 %w/w

**Sample size** 0.3g

**Summary of analytical procedure** Fill the vessel with 300 mL of solvent mixture. Calibrate with 25uL of DI water. Start titration in Bgrd-Samp-Bgrd mode. Accurately weigh 1 capsule and enter total weight when display requests. Add starch from the capsule to the titration vessel. Reweigh the empty gelatin capsule and enter as holder weight when display requests.

### Method Parameters

**Instrument** Turbo 2 Volumetric Karl Fischer

**Volumetric Reagent** HYDRANAL®-Composite 5 RDH 34805\*

**Volumetric Solvent** 1:1 Methanol:Formamide (Karl Fischer grade)

**Endpt/Step Level** 8/ +5

**Endpt Time** 10 sec.

**Stirrer/blender speed** 3

**Extract time** 10 min.\*\*

**Mode of analysis** Bgrd-Samp-Bgrd

**Vaporizer oven** N/A

**Result units** % w/w

### Results

**# of Trials** 1

**Mean Value** 8.5 %

**C.V.** N/A

**Comments** \*Or use AquaStar\* Composite 5 from Fisher (MAX1698A-6). \*\*Turbo time can be reduced to 5 minutes if the titration is carried out at 50°C. Use the Turbo2 Water-Jacketed Vessel with circulating bath to perform calibration and analysis at 50°C.

**Cross-reference log #** KF1142