

Copper Base

ARL 9900 Series with IntelliPower™
Simultaneous-Sequential XRF Spectrometer


Key Words

- ARL 9900 - 3600W
- Brass
- Bronze
- Copper
- XRF

Copper alloys

Pure copper is mixed with other elements to produce a wide range of alloys:

- Brasses (Cu-Zn)
- Bronzes (Cu-Sn)
- Gun metal (Cu-Sn-Zn) used for casting due to its excellent fluidity, for valves, taps and water fittings
- Manganese bronzes (Cu-Zn-Mn) are not true bronzes as Sn is not the essential constituent. Used for rudders, propellers and ship fittings
- Cupro-nickel (Cu-Ni) used for coins, tubes, wires, electrical resistances and thermocouples
- Nickel silver (Cu-Ni-Zn) used for marine applications, car radiators and fittings



Brasses

One of the most important non-ferrous engineering alloys. They cover a large range of physical properties and their applications are multiple.

Copper is capable of holding about 39 % of zinc in solid solution. Alloys containing less than 39 % Zn are known as α -brasses. This brass (70/30) is very widely employed for cartridge, cases, condenser tubes, etc.

From 39 % to 46 % of Zn, a β -solid solution gives $\alpha\beta$ -brass. This alloy (60/40) is found in extensive engineering applications with enhanced corrosion resistance used in marine applications. Alloys containing more than 49 % Zn are very hard and only used for brazing brasses.

Bronzes

Bronzes are used for bearing and gears. There are a wide variety of bronze alloys containing various elements such as phosphorus, beryllium, silicon, etc.:

- Leaded bronzes (Cu-Sn-Pb) are mainly used for bearings due to their good wear resistance
- Aluminum bronzes (Cu-Al) have good corrosion resistance but are difficult to cast. Used for pumps rods, diecastings, etc.
- Silicon bronzes (Cu-Sn-Ni) have high electrical conductivity and are used for wires
- Phosphor bronzes (Cu-Sn-P) are used for valves, bearings and gears. Cold worked into rods and sheets

Typical performance

ELEMENT		CONCENTRATION %	PRECISION (1 SIGMA) %	LOD (3 SIGMA) PPM
COUNTING TIME			30 s.	30 s.
Al	Aluminum	0.01	0.0006	6
		0.1	0.004	
		5	0.007	
		10	0.01	
Ag	Silver on gonio + PBF	0.01	0.001	30
		0.05	0.0015	
As	Arsenic K β	0.01	0.0003	9
		0.06	0.0004	
Bi	Bismuth on gonio	0.01	0.0004	21
		0.04	0.0007	
Cd	Cadmium on gonio + PBF	0.01	0.0017	45
		0.025	0.002	
Cr	Chromium	0.01	0.0001	3
Co	Cobalt	0.005	0.00003	1
Cu	Copper	50	0.01	--
		65	0.015	
		80	0.02	
Fe	Iron	0.03	0.00005	2
		0.1	0.0002	
		2	0.0005	
		5	0.002	
Mn	Manganese	0.01	0.0001	2
		1	0.0005	
		5	0.001	
Ni	Nickel	0.01	0.0001	4
		0.2	0.0003	
		1	0.001	
		5	0.003	
		10	0.004	
		30	0.006	
P	Phosphorus Ge 111	0.01	0.0003	17
		1	0.0004	
Pb	Lead L β	0.01	0.0002	4
		1	0.001	
		5	0.005	
S	Sulphur	0.02	0.0004	4
		0.1	0.001	
Sb	Antimony	0.01	0.0002	7
		0.1	0.0006	

(cont. on back page)

ELEMENT		CONCENTRATION %	PRECISION (1 SIGMA) %	LOD (3 SIGMA) PPM
COUNTING TIME			30 s.	30 s.
Se	Selenium on gonio	0.01	0.0004	7
Si	Silicon on gonio	0.01	0.0001	8
		0.05	0.0004	
		0.1	0.0007	
Sn	Tin	0.02	0.0002	5
		0.5	0.0003	
		5	0.0014	
		10	0.002	
Te	Tellurium on gonio	0.01	0.0012	30
		0.03	0.0015	
Zn	Zinc	0.2	0.0005	8
		1	0.0007	
		5	0.0015	
		20	0.008	
		40	0.011	

This list will be updated as improvements are announced. Please contact your nearest Thermo office for the most recent values.

The values listed in the general table are calculated from the following formulae:

$$\text{Precision (1 sigma)} = \sqrt{\frac{P + B}{Q T}}$$

and

$$\text{LOD} = \text{limit of detection (3 sigma)} = 3 \sqrt{\frac{B}{Q T}}$$

where

- Q = counts per second per 1 % element
- B = background equivalent concentration
- T = time of analysis in seconds
- P = the percentage concentration at which the sigma value is calculated.

The limits of detection are determined with low alloy copper samples. In practice, these values should agree with values measured practically using the mathematical derivation of standard deviation (1 sigma).

$$1 \text{ sigma} = \pm \sqrt{\frac{\sum(X - \bar{X})^2}{n - 1}}$$

where

- \bar{X} = the individual readings
- \bar{X} = the arithmetic mean of the individual values
- n = number of determination (normally ≥ 10)

For guaranteed values, the values of precision and LODs should be multiplied by factor 1.5.

Instrument calibration

It should be stressed that an XRF spectrometer is a very accurate comparator, but the accuracy of the final analysis is entirely dependent on the quality of the standards used for calibration. The Thermo Scientific ARL 9900 XRF spectrometer, can be factory calibrated for various copper alloys or a general copper program using commercially available standards. It can also be calibrated on-site using well analyzed samples from the user.



Conclusion

The long experience of our company in metals analysis comes from an installed base of over 1'000 XRF spectrometers world-wide. The ARL 9900 is the answer to your metallurgical analysis needs; whether they be incoming material control, metal QC or production analysis. Working 24 hours a day and 7 days a week, the ARL 9900 delivers dependable performance year after year. The high performance of the ARL 9900 will meet your analytical needs today and in the future.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa
+43 1 333 5034 127

Australia
+61 2 8844 9500

Austria
+43 1 333 50340

Belgium
+32 2 482 30 30

Canada
+1 800 530 8447

China
+86 10 5850 3588

Denmark
+45 70 23 62 60

Europe-Other
+43 1 333 5034 127

France
+33 1 60 92 48 00

Germany
+49 6103 408 1014

India
+91 22 6742 9434

Italy
+39 02 950 591

Japan
+81 45 453 9100

Latin America
+1 608 276 5659

Middle East
+43 1 333 5034 127

Netherlands
+31 76 587 98 88

South Africa
+27 11 570 1840

Spain
+34 914 845 965

Sweden/Norway/Finland
+46 8 556 468 00

Switzerland
+41 21 694 71 11

UK
+44 1442 233555

USA
+1 800 532 4752

www.thermo.com



Thermo Electron SA, Ecublens, Switzerland is ISO certified.

AN41401_E 01/08C

©2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.