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- XI. Referenzen**
- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
  - Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
  - NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS Document H21-A4. NCCLS, Wayne, PA, 2003.
  - Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
  - Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
  - Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
  - Palaereti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise.* Thromb Haemostasis 58:905-910, 1987.
  - Stabilitätsdaten in DHF.
  2002. Daten in 510 (K)-Datei.

Bestellinformationen		
<b>Kat.- Nr.</b>	<b>Beschreibung</b>	<b>Inhalt</b>
100354	Thromboplastin-DS	10 x 4 mL
100362	Thromboplastin-DS	10 x 10 mL
100353	Thromboplastin-DS	10 x 2 mL

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- Technique:**
- pH will increase if plasma is open to air. Store samples stoppered in plastic or siliconized glass.
  - Plasma held at 4-8°C may undergo cold activation resulting in significant shortening of the PT.<sup>4</sup>
  - Thromboplastin-DS was designed to work at 37°C ± 0.5°C. Frequently check the temperature of all heating elements.
  - All labware must be clean and free of trace amounts of detergents.
  - Always follow instrument manufacturer’s instructions for proper maintenance.
- Interfering Substances:**
- Sodium oxalate, EDTA, and heparin are not suitable anticoagulants.
  - The PT may be prolonged by substances such as oral contraceptives, corticosteroids, EDTA, asparaginase, clofibrate, erythromycin, ethanol, tetracycline, and anticoagulants such as heparin and warfarin.<sup>5</sup>
  - The PT may be shortened by substances including antihistamines, butabarbital, caffeine, oral contraceptives, phenobarbital, and vitamin K.<sup>5</sup>

**IX. Expected Values**

In multi-center studies, when Thromboplastin-DS was evaluated on a normal population, the following results were obtained:<sup>9</sup>

Instrument	PT Mean(sec)	Range (+/-2SD)	N
MLA™ Electra 1000CT™	13,2	11,4–15,0	40
MLA™ Electra 900CT™	13,7	12,4–15,0	20
IL ACL™ 300/3000+	10,5	8,9–12,1	61
Amelung KC 10™	12,7	9,3–14,2	20
Pacific Hemostasis ThromboScreen 400C	13,5	12,2–14,8	38
Pacific Hemostasis ThromboScreen 200	13,5	12,0–15,1	60

These values are intended as a guideline only. Each laboratory should establish a Normal Reference Range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory. The NRR should be reestablished or at least verified when changing lot numbers of the same reagent.<sup>3,6</sup> A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.

The clotting time of abnormal plasmas will depend on the ISI of the reagent lot in use.

**X. Performance Characteristics**

**Precision:** Precision of Prothrombin Time results is dependent on many factors, such as the instrument, technique and the reagent in use. Thromboplastin-DS precision was assessed by testing a normal and abnormal plasma on several different instruments. A summary of the results follows.<sup>10</sup>

**Summary of Within-run Precision Studies, % CV (N = 20)**

Sample	MLA Electra 1000C	ThromboScreen 400C	ThromboScreen 200	Amelung KC10
<b>Normal</b>	1,1 <span> </span> %	1,9 <span> </span> %	1,9 <span> </span> %	2,9 <span> </span> %
<b>Abnormal</b>	2,8 <span> </span> %	2,5 <span> </span> %	2,3 <span> </span> %	1,1 <span> </span> %

**Sensitivity:** Thromboplastin-DS detects deficiencies in the extrinsic pathway as determined by the Prothrombin Time test. Factor sensitivity testing was performed by diluting pool normal plasma with factor deficient plasmas such that the final factor concentration ranged from 10-100%. PT testing of these samples was performed on the MLA-1000C instrument.<sup>11</sup>

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**Correlation:** Correlation studies were performed against two other sensitive thromboplastin reagents by performing PT testing on normal and abnormal samples. Testing was performed on the Stago STA instrument.<sup>12</sup>

	PT Correlation	INR Correlation
<b>Thromboplastin DS vs. Reagent A, N = 49</b>	R = 0,98	R = 0,98
<b>Thromboplastin DS vs. Reagent B, N = 49</b>	y = 1,16x + 1,30	y = 0,89x + 0,05
	R = 0,95	R = 0,95
	y = 1,01x + 2,20	y = 0,82x + 0,10

**XI. Referenzen**

- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic.* Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range.* Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline.* NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time.* Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests.* Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy.* Arch Inter Med. 146:462-472, 1986.
- Palaereti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise.* Thromb Haemostasis 58:905-910, 1987.
- Stability data found in DHF.
2002. Data found in 510(K) file.

Ordering Information		
<b>Cat. No.</b>	<b>Description</b>	<b>Contents</b>
100354	Thromboplastin-DS	10 x 4 mL
100362	Thromboplastin-DS	10 x 10 mL
100353	Thromboplastin-DS	10 x 2 mL

**FISHER DIAGNOSTICS® LIMITED WARRANTY**
Fisher Diagnostics (FD) warrants to the purchaser only that FD products will perform as described on their labeling and product literature. Purchaser must determine the suitability of FD products for their specific applications. FD’s sole obligation will be, at its option, to either replace a non-conforming or defective product, or return the purchase price. FD DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE. Neither FD nor its affiliates shall, in any event, be liable for incidental or consequential loss or damage.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. ACL is a registered trademark of Instrumentation Laboratory. MLA Electra 900C and MLA Electra 1000C are registered trademarks of Instrumentation Laboratory. Amelung KC10 is a registered trademark of Trinity Biotech

## Pacific Hemostasis® Thromboplastin-DS

- I. Intended Use**

Pacific Hemostasis Thromboplastin-DS is intended for use in performing the one-stage prothrombin time (PT) test and PT-based factor assays.

**II. Summary and Principles**

The PT is used as a screening tool and as a quantitative test for coagulation factors in the extrinsic and common pathways. This test will be prolonged in patients with acquired or congenital disorders that reduce the activity of factors I (fibrinogen), II (Prothrombin), V, VII, and X. The PT is also widely used to monitor oral anticoagulant therapy.<sup>1,2</sup> Oral anticoagulants reduce the activity of vitamin-K dependent clotting factors (II, VII, IX, X, Protein C, and Protein S), and the PT is prolonged as a result.

The one-stage PT measures the clotting time of plasma after adding a source of tissue factor (thromboplastin) and calcium. The recalcification of plasma in the presence of tissue factor generates activated Factor Xa (F:Xa). F:Xa in turn activates Prothrombin to thrombin, which converts fibrinogen to an insoluble fibrin clot.

**III. Reagent**

For *in vitro* diagnostic use.

**Composition:** <2 % rabbit brain tissue, 0,013 % sodium azide, 5% buffers, salts, and stabilizers
Store unopened vials at 2-8°C. Reconstitute with preservative-free distilled/deionized water according to the Thromboplastin-DS vial label, swirl gently and let the vial stand undisturbed for 15 minutes at room temperature. Do not invert the vial or mix vigorously. Reconstitution fluid is available if the quality of water is questionable. After reconstitution, the reagent when stored stoppered is stable for 7 days at 2-8°C, 8 hours at 37°C, and for 24 hours at room temperature (15-25°C).<sup>8</sup> Store at 2-8°C when not in use. **Do not freeze.**

Mix gently before each use. Provide some mechanism, such as a magnetic stirrer, to maintain adequate suspension during use. Lack of vacuum in vials, erratic results, quality control values outside established ranges, or product color variations could indicate deterioration. However, poor performance could also be due to other factors within the test system.

**Warning:** Thromboplastin-DS contains sodium azide. Sodium azide under acidic conditions yields hydrazoic acid, an extremely toxic compound. Azide compounds should be diluted with running water before being discarded. Upon disposal, azide compounds should be flushed with large volumes of water. These precautions are recommended to avoid deposits in metal pipes in which explosive conditions may develop.

**IV. Specimen Collection**

3.2% (0.109M) trisodium citrate anticoagulant is recommended for coagulation assays. Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90% of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22-24°C. If testing is not completed within 24 hours, plasma should be frozen at –20°C for up to two weeks or –70°C for up to 6 months. For more details on specimen collection and storage, see NCCLS Document H21-A4.<sup>3</sup>

- Do not delay mixing the blood with anticoagulant.
- Avoid foaming the specimen.
- Use only plastic or siliconized glass containers.
- Turbid, icteric, lipemic, or hemolyzed specimens may generate erroneous results.
- Freezing and thawing plasma that contains residual cells will generate damaged cell membranes that can affect results.
- Acute inflammatory reactions can shorten PT results because of elevated fibrinogen.
- Plasma samples with hematocrits outside the range of 20-55% may be improperly anticoagulated and should be adjusted appropriately.

**V. Test Procedure**

**Materials Provided:** Thromboplastin-DS Reagent, 10 x 4 mL, 10 x10 mL, or 10 x 2 mL

**Materials Required, But Not Provided:**

Distilled or Deionized water, or Pacific Hemostasis Reconstitution Fluid

Stopwatch or timer

Precision pipette: 0.1 and 0.2 mL

Normal and abnormal controls such as Pacific Hemostasis Coagulation Control Plasmas, Level 1, 2, and 3

Thromboplastin-DS is suitable for use with manual, mechanical, photo-optical, nephelometric or other means of clot detection. Follow manufacturer’s recommendations for proper use of instrumentation. For manual assays:

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**XI. Bibliografía**

- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic*. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range*. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline*. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time*. Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests*. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy*. Arch Inter Med. 146:462-472, 1986.
- Palaereti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise*. Thromb Haemostasis 58:905-910, 1987.
- Los datos de la estabilidad proceden de DHF.
- 12- Información procedente de un archivo de 510(K).

	Información para pedidos		
Nº de ref.	Descripción	Contenido	
100354	Thromboplastin-DS	10 x 4 mL	
100362	Thromboplastin-DS	10 x 10 mL	
100353	Thromboplastin-DS	10 x 2 mL	

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All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

ACL is a registered trademark of Instrumentation Laboratory.
MLA Electra 900C and MLA Electra 1000C are registered trademarks of Instrumentation Laboratory.
Amelung KC10 is a registered trademark of Trinity Biotech

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- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral traitement anticoagulant: Experience with an Anticoagulation Clinic*. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range*. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline*. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time*. Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests*. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy*. Arch Inter Med. 146:462-472, 1986.
- Palaereti, G., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise*. Thromb Haemostasis 58:905-910, 1987.
- Données de stabilité provenant du DHF.
- 12- Données figurant dans le dossier 510(k).

	Pour commander		
Réf. Cat.	Description	Contenu	
100354	Thromboplastin-DS	10 x 4 mL	
100362	Thromboplastin-DS	10 x 10 mL	
100353	Thromboplastin-DS	10 x 2 mL	

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All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.
ACL is a registered trademark of Instrumentation Laboratory.
MLA Electra 900C and MLA Electra 1000C are registered trademarks of Instrumentation Laboratory.
Amelung KC10 is a registered trademark of Trinity Biotech

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#### XI. Bibliografia

- Erichetti, A.M., Holden, A., Ansell, J.: *Management of Oral Anticoagulant Therapy: Experience with an Anticoagulation Clinic*. Arch Inter Med. 144:1966-68, 1984.
- Hirsh, J., Dalen, J.E., Deykin, D., Poller, L.: *Oral Anticoagulants: Mechanisms of Action, Clinical Effectiveness, and Optimal Therapeutic Range*. Chest 102(Suppl): 312S-316S, 1992.
- NCCLS: *Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline*. NCCLS document H21-A4. NCCLS, Wayne, PA, 2003.
- Palmer, R.N., Gralnick, H.R.: *Inhibition of the Cold Activation of Factor VII and the Prothrombin Time*. Am J Clin Path. 81: 618-622, 1984.
- Young, D.S., Thomas, D.W., Friedman, R.B., et al: *Effect of Drugs on Clinical Laboratory Tests*. Clin Chem 18:1041, 1972.
- Dalen, J.E., Hirsh, J.: *American College of Chest Physicians and the National Heart, Lung, and Blood Institute National Conference on Antithrombotic Therapy*. Arch Inter Med. 146:462-472, 1986.
- Palaueri, C., Coccheri, S., Poggi, M., et al: *Oral Anticoagulant Therapy Control: Evidence that the INR Expression Improves the Interlaboratory Comparability of Results. The Bologna Oral Anticoagulant Control Exercise*. Thromb Haemostasis 58:905-910, 1987.
- I dati relativi alla stabilità sono stati ottenuti da DHF.
- 9-12. Informazioni ottenute dal file 510(k).

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Fisher Diagnostics  
a division of Fisher Scientific Company, LLC  
a part of Thermo Fisher Scientific Inc.  
Middletown, VA 22645-1905 USA  
Phone: (800) 528-0494  
(540) 869-3200  
Fax: (540) 869-8132

840587

Rev. R1

**Thermo**  
S C I E N T I F I C

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MDCI Ltd.