

Cellomics[®] Whole Cell Stains Blue, Green, Orange and Red

For High Content Screening and Fluorescence Microscopy

1998.1

| Number | Description | Storage | Ex/Em |
|---------|---------------------------------|---------|------------|
| 8403501 | Whole Cell Stain Blue, 1 × 96 | 4°C | 350/450 nm |
| 8403502 | Whole Cell Stain Blue, 5 × 96 | | |
| 8403201 | Whole Cell Stain Green, 1 × 96 | -20°C | 493/518 nm |
| 8403202 | Whole Cell Stain Green, 5 × 96 | | |
| 8403301 | Whole Cell Stain Orange, 1 × 96 | 4°C | 550/568 nm |
| 8403302 | Whole Cell Stain Orange, 5 × 96 | | |
| 8403401 | Whole Cell Stain Red, 1 × 96 | 4°C | 654/673 nm |
| 8403402 | Whole Cell Stain Red, 5 × 96 | | |

Storage: Store the product as indicated above. Keep stains away from light and moisture. See the **Solution Preparation** section for storage and stability of prepared solutions.

Warning: Please completely read these instructions and the accompanying material safety data sheets before using this product. Cellomics Reagents are not for diagnostic use in humans or animals.

Introduction

The Cellomics Whole Cell Stains provide excellent staining for HCS assays and fluorescence microscopy. These stains are intense, highly photostable and match the output wavelengths of common fluorescence instrumentation (Figure 1). Effective image analysis in HCS cell-based assays and fluorescence microscopy requires fluorescent labeling of the entire cell. In these assays, the cellular primary object is used to identify and count individual cells and define the cell region in which the image analysis is applied. The primary object might be a major cellular component, such as the nucleus, a large organelle or the whole cell. When the whole cell is the primary object, high-quality Cellomics Whole Cell Stains effectively distinguish intact cells from bordering cells (Figure 2).




Figure 1. Excitation and emission spectra for the Cellomics Whole Cell Stains: Panel A blue stain in methanol, Panel B green stain in PBS, Panel C orange stain in PBS and Panel D red stain in PBS.




Figure 2. Image of human U2OS cells stained with Whole Cell Stain Green (with 20X objective). U2OS cells were stained for 30 minutes as described in the protocol. Whole Cell Stain Blue, Orange and Red stain the cell correspondingly.

Additional Materials Required

- Phosphate-buffered saline (PBS) (Thermo Fisher Scientific, Pierce Product No. 28372)
- Dimethylsulfoxide (DMSO)

Solution Preparation

| | |
|---|--|
| Whole Cell Stain Stock Solution | For Whole Cell Stain Green, Orange and Red, add 50 μ l of DMSO to the 1 \times 96 vial or 250 μ l DMSO to the 5 \times 96 vial and mix. If there is any stain lodged under the cap, reconstitute it carefully using a pipette. For Whole Cell Stain Blue, add 100 μ l of DMSO to the 1 \times 96 vial or 500 μ l DMSO for the 5 \times 96 vial and mix. Discard unused Whole Cell Stain Stock Solution or for short-term storage (< 1 week), store at -20°C. |
| Whole Cell Stain Working Solution (per 96-well plate) | For Whole Cell Stain Green, Orange and Red, dilute 50 μ l of the Whole Cell Stain Stock Solution to 5 ml with PBS and mix well. For Whole Cell Stain Blue, dilute 100 μ l of the stain stock solution to 5 ml with PBS and mix well. Note: For best results, use the whole cell stain working solution immediately; however, it remains effective for at least 3 hours at room temperature. After staining, discard any unused stain working solution. |

Whole Cell Staining Procedure

- Do not allow the wells to dry at any time during the protocol.
 - Allow all vials to warm to room temperature before opening.
 - Perform all steps at room temperature unless indicated otherwise.
 - Using conditions other than those specified in these instructions may require protocol optimization.
1. Fix cells with formaldehyde at room temperature. If desired, probe with specific antibodies or dyes before using the Whole Cell Stain.
Note: Cell permeabilization (e.g., 0.1% Triton X-100 in PBS for 15 minutes) may enhance staining intensity.
 2. Rinse cells 2-3 times with PBS.
 3. Add 50 μ l of the Whole Cell Stain Working Solution to each well.
 4. Incubate cells in the dark for 30 minutes at room temperature. Optimal staining time varies from 10 minutes to 3 hours depending on cell type.
 5. Aspirate the Whole Cell Stain Working Solution and wash three times with PBS.
 6. Seal plate and evaluate on a HCS reader. For fluorescence microscopy, observe the stained cells using the proper excitation and emission filters.
 7. Store sealed plates in the dark at 4°C.

Recommendations for Automation

- **Plating Cells:** To improve the uniformity and throughput of plating cells, use a liquid-handling system such as a Titertek[®] Multidrop Dispenser.
- **Dead Volumes:** Every piece of automation instrumentation has a non-recoverable dead volume associated with it. Be aware of these dead volumes, priming volumes and rinsing volumes when calculating your reagent requirements.
- **Nonspecific Binding:** Because of the potential of reagent interaction with large surface areas inherent to tubing, syringes and peristaltic pumps, pre-priming with reagents or pre-coating with protein blockers may be warranted.
- **Mixing:** Gentle mixing may be required when adding a DMSO-based solution to keep overly concentrated solutions from lying on top of the cell layer. Be careful not to dislodge cells or beads during mixing procedures.
- **Cell Washing:** Use an automated plate washer designed to gently wash attached cells. Be careful not to dislodge cells or beads during cell washing.
- **Incubation:** Minimize the time when plates with live cells are out of a controlled CO₂ environment. For best results, use an automated incubator to deliver plates to a pipetting deck.
- **Exposure:** Minimize operator exposure to fixative by some form of containment. Some reagents and compounds are light-sensitive; be aware of these constraints when scaling up for an automated run.

The listed BioApplications Software Modules are protected by U.S. patent #5,989,835 and other patent pending. Titertek[®] is a registered trademark of Titertek Instruments, Inc.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current versions of product instructions are available at www.thermo.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2007 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.