EN

**CRP Control**

REF 981251 5 x 1 ml

**THIS PACKAGE INSERT IS APPLICABLE FOR USE OUTSIDE THE US.****INTENDED USE**

CRP Control is intended for use on Konelab and Indiko Clinical Chemistry Analyzers as a quality control serum to monitor trueness of Thermo Fisher Scientific Oy defined methods CRP (codes 981699 and 981933) and CRP Plus (codes 981794 and 981934).

Any reference to the Konelab systems also refers to the T Series.

**COMPOSITION**

CRP Control serum is a liquid human based reference preparation for immunoturbidimetric methods. This product contains < 0.1% sodium azide as a preservative.

**PRECAUTIONS**

For *in vitro* diagnostic use only. This product has been tested, and found negative for HBsAg and for antibodies to HIV 1/2- and HCV. However, the product must be handled just as carefully as the patient samples.

**HANDLING**

Liquid, ready for use. Allow the vial to attain the room temperature prior to opening. Agitate gently avoiding foam formation.

It is recommended to let the control be max. for two hours on the analyzer.

**STORAGE AND STABILITY**

CRP Control is stable as follows:

|                     |   |
|---------------------|---|
| In an unopened vial | Until the expiration date when stored at 2...8 °C |
| After opening       | 2 months at 2...8 °C                              |

Store the control tightly capped when not in use.

**LIMITATIONS**

Do not use any vial that has visible evidence of microbial growth.

**ASSIGNED VALUES AND TRACEABILITY**

The concentration of CRP is lot-specific. The values are listed in the analyzer specific Value Sheets. Assigned values for CRP is given for Konelab and Indiko defined methods and analyzers.

It is recommended that each laboratory establish its own means and acceptable ranges. The control intervals and limits must be adapted to the individual laboratory requirements. Primary reference: ERM-DA470.

**CHANGING CONTROL VALUES**

Check always the lot number of the control in your hand. Note that an additional letter, e.g. A, at the end of numeric lot code specifies the packaging time.

**Konelab**

If the lot number is new, change the respective values in the parameter list manually.

**Indiko**

If the lot number is new, you can read the respective values either from a barcode or from a file, or you can enter the values manually. For more information, see the *Operation Manual* for Indiko.

**MANUFACTURED BY**

Thermo Fisher Scientific Oy  
Clinical Diagnostics Finland  
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland  
Tel. +358 9 329 100, Fax +358 9 3291 0300  
www.thermoscientific.com

**Date of revision (yyyy-mm-dd)**

2011-08-29

**Changes from previous version**

The line beside the text indicates changes.

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