EN

**NORTROL**

REF 981043 10 x 5 ml

**THIS PACKAGE INSERT IS APPLICABLE FOR USE OUTSIDE THE US. IN THE US, THIS PACKAGE INSERT IS APPLICABLE FOR KONELAB 20XT(i), KONELAB PRIME 60(i), T60.**

**INTENDED USE**

Nortrol is a control serum to monitor trueness and precision of the analytes listed in the separate Nortrol value sheet. The given values are valid for Konelab and Indiko Clinical Chemistry Analyzers using methods defined by Thermo Fisher Scientific Oy. Any reference to the Konelab systems also refers to the T Series.

**PRODUCT INFORMATION**

Kit code 981043 contains 10 vials of 5 ml Nortrol control.

**COMPOSITION**

Nortrol is prepared from human serum. This product contains human and non-human enzymes as well as non-protein constituents. Bacteriostatic agents have been added.

**PRECAUTIONS**

For *in vitro* diagnostic use only. Exercise the normal precautions required for handling all laboratory reagents. This product has been tested by FDA approved methods, and found negative for HBsAg and for antibodies to HIV 1/2- and HCV. However, the product must be handled just as carefully as the patient samples.

**HANDLING**

- To reconstitute Nortrol,
- 1) Remove the cap and carefully free the rubber stopper avoiding the loss of lyophilizate.
  - 2) Add exactly 5.00 ml of distilled or \*deionized water.
  - 3) Close the vial and allow to stand for 30 min at room temperature.
  - 4) Dissolve the contents completely by swirling gently. Do not shake.
- \* For electrolyte measurements, preferably use deionized water.

It is recommended to let the control be max. for two hours on the analyzer.

**STORAGE AND STABILITY**

Nortrol is stable as follows:

In an unopened vial	Until the expiry date when stored at 2...8°C
After reconstitution	7 days at 2...8°C, or 30 days at -20°C

**Exceptions (after reconstitution)**

Acid Phosphatase	Must be analyzed immediately; do not freeze
Alkaline phosphatase	3 days at 2...8°C

Store the controls tightly capped when not in use. Prolonged exposure to light may alter the bilirubin concentration.

**LIMITATIONS**

Do not use any vial that has visible evidence of microbial growth.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

**ASSIGNED VALUES AND TRACEABILITY**

The concentrations of the components are lot-specific. The values are listed in the analyzer specific Value Sheets.

Assigned values and traceability for the components of Nortrol are given for the methods defined by Thermo Fisher Scientific Oy. These values for trueness and precision controls are traceable to the reference materials listed in a separate traceability certificate, which is provided on request.

It is recommended that each laboratory establish its own means and acceptable ranges. The control intervals and limits must be adapted to the individual laboratory requirements. Always follow the local, state and federal regulations in performing quality control (QC).

**CHANGING CONTROL VALUES**

Check always the lot number of the control in your hand. Note that an additional letter, e.g. A, at the end of numeric lot code specifies the packaging time.

**Konelab**

If the lot number is new, change the respective values in the parameter list manually.

**Indiko**

If the lot number is new, you can read the respective values either from a barcode or from a file, or you can enter the values manually. For more information, see the *Operation Manual* for Indiko.

**MANUFACTURED BY**

Thermo Fisher Scientific Oy  
Clinical Diagnostics Finland  
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland  
Tel.: +358 9 329 100, Fax: +358 9 3291 0300  
www.thermoscientific.com

**Date of revision (yyyy-mm-dd)**

2011-08-10

**Changes from previous version**

The line beside the text indicates changes.

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