EN

**Lipotrol HDL/LDL
Abnormal**

REF 981907 5 x 3 ml

**THIS PACKAGE INSERT IS APPLICABLE
FOR USE OUTSIDE THE US.****INTENDED USE**

Lipotrol HDL/LDL Abnormal is a control serum to monitor trueness and precision of HDL-Cholesterol and LDL-Cholesterol methods defined by Thermo Fisher Scientific Oy on Konelab and Indiko Clinical Chemistry Analyzers. Any reference to the Konelab systems also refers to the T Series.

COMPOSITION

5 vials, each with lyophilizate for 3 ml of control.
Lipotrol HDL/LDL Abnormal is a human based lyophilized control serum.

PRECAUTIONS

For *in vitro* diagnostic use only. This product has been tested, and found negative for HBsAg and for antibodies to HIV 1/2- and HCV. However, the product must be handled just as carefully as the patient samples.

HANDLING

To reconstitute Lipotrol HDL/LDL Abnormal,

- 1) Remove the cap and carefully free the stopper.
- 2) Add exactly 3.00 ml of distilled or deionized water.
- 3) Close the vial and allow to stand for 30 minutes.
- 4) Dissolve the contents by swirling gently, avoiding foam formation.

It is recommended to let the control be max. for two hours on the analyzer.

STORAGE AND STABILITY

Lipotrol HDL/LDL Abnormal is stable as follows:

In an unopened vial	Until the expiration date when stored at 2...8 °C
After reconstitution	5 days at 2...8 °C, or 4 weeks at -20 °C when frozen once

Store the control tightly capped when not in use.

LIMITATIONS

Do not use any vial that has visible evidence of microbial growth.

ASSIGNED VALUES

The concentrations of HDL- and LDL-Cholesterol are lot-specific. The values are listed in the analyzer specific Value Sheets.

Assigned values for the components of Lipotrol HDL/LDL Abnormal are given for the methods defined by Thermo Fisher Scientific Oy.

CHANGING CONTROL VALUES

Check always the lot number of the control in your hand. Note that an additional letter, e.g. A, at the end of numeric lot code specifies the packaging time.

Konelab

If the lot number is new, change the respective values in the parameter list manually.

Indiko

If the lot number is new, you can read the respective values either from a barcode or from a file, or you can enter the values manually. For more information, see the *Operation Manual* for Indiko.

MANUFACTURED BY

Thermo Fisher Scientific Oy
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland
Tel. +358 9 329 100, Fax +358 9 3291 0300
www.thermoscientific.com

Date of revision (yyyy-mm-dd)

2011-08-26

Changes from previous version

The line beside the text indicates changes.

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