

Concentration of Retrovirus Using Thermo Scientific Sorvall® WX Ultra Centrifuge and Thermo Scientific SureSpin™ 630 Swinging Bucket Titanium Rotor

Stephanie R. Noles, PhD, Thermo Fisher Scientific, Laboratory Equipment Division, Associate Product Manager

Key Words

- Ultracentrifugation
- Retrovirus Concentration
- Swinging Bucket Rotor



Introduction

Protein function can be determined by using gain- or loss- of function experiments including the induction of exogenous gene expression in multiple cell lines and *in vivo* models. The utilization of viral transduction to induce exogenous gene expression is becoming an increasingly common technology. This technology is more efficient than transfection as you can selectively target cell populations by using multiple virus types (ie, retrovirus, adenovirus, lentivirus), gene expression is increased in difficult to transfect cell types, and gene expres-

sion in infected cells is permanent (Cepko and Pear, 1996). Multiple methods can be used to conduct viral transduction including incubation of target cell populations with dilute or concentrated viral supernatant. Furthermore, concentrated virus can be used for the rapid induction of genetic expression *in vivo* which can eliminate the need for costly and time-consuming transgenic mouse models (Gaiano et al., 1999). The concentration of virus can be completed by ultracentrifugation of viral supernatant harvested from transfected packaging cell lines.

The procedure described below outlines the process for the production and concentration of retrovirus using Thermo Scientific Sorvall WX Ultra Centrifuge and Thermo Scientific SureSpin 630 Swinging Bucket Titanium Rotor. By incubating cortical progenitors with concentrated retrovirus produced by this method, we were able to selectively target gene expression (identified by green fluorescent protein, GFP) to dividing neural progenitors, a primary cell line that has difficulty with exogenous gene expression.

Procedure

Production of dilute virus

- Transfect GP2-293 derived packaging cell line (Clontech) at 90% confluency with 5µg each of pCLE-IRES2-eGFP (Anjen Chenn, Northwestern University, Chicago, IL) and pVSV-G (Clontech) using Lipofectamine™ and PLUS™ Reagent (Invitrogen)
- Approximately 18hrs post-transfection, replace the culture media with fresh media according to manufacturer's instructions
- Place transfected GP2-293 cells at 32°C and 5% CO₂ for virus production
- Harvest virus by removing viral supernatant from confluent monolayer of transfected cells 5 to 6 times at 12 hr intervals and placing into conical tube
- Centrifuge conical tube in Thermo Scientific lowspeed centrifuge at 200 x g for 5 mins at room temperature
- Syringe tube filter viral supernatant through a 45µm filter into new conical
- Freeze each filtered harvest at -80°C




FIGURE 1. Retrovirally infected neural progenitor cells identified by GFP expression. Blue, nuclei, Green, Green Fluorescent Protein.

Concentration of dilute virus

- Thaw viral harvests at room temperature
- Remove viral harvest(s) into PA Thin-Walled Tubes (36 mL, cat #: 03141)
- Place balanced tubes into Thermo Scientific SureSpin 630 (36 mL) Swinging Bucket Titanium Rotor
- Centrifuge tubes in WX Ultra Centrifuge at 112,700 x g (25,000rpm) at 4°C for 1.5 hrs
- Remove supernatant from tubes, leaving viral pellet at the bottom
- Resuspend viral pellet in 40 µL sterile PBS overnight at 4°C and aliquot into small volume working stocks in microcentrifuge tubes
- Store aliquots at -80°C

Discussion

The procedure above describes a method to produce concentrated retrovirus using the WX Ultra Centrifuge and SureSpin 630 Swinging Bucket Rotor. Using concentrated retrovirus we were able to selectively infect dividing neural progenitor cells as detected by GFP expression in small volume without disturbing the culture conditions of the cells (Figure 1, Noles and Chenn, 2007). The above procedure can also be used to produce concentrated retrovirus to induce rapid exogenous gene expression in multiple *in vivo* models (ie, mouse and chick) without the need to create expensive and time consuming transgenic mouse models.

References

- Cepko, Constance and Warren Pear. "Transduction of Genes Using Retrovirus Vectors." *Current Protocols in Molecular Biology* (1996) 9.9.1-9.9.16
- Gaiano N, Kohtz JD, Turnbull DH, and Fishell, G. (1999). A method for rapid gain-of-function studies in the mouse embryonic nervous system. *Nat Neurosci.* (9):812-9.
- Noles, SR and Anjen Chenn (2007). Cadherin inhibition of beta-catenin signaling regulates the proliferation and differentiation of neural precursor cells. *Mol. Cell. Neurosci.* (35): 549-58.



Thermo Scientific
SureSpin 630
Swinging Bucket
Titanium Rotor

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

North America:
USA / Canada
+1 800 553 0039

Europe:
Austria
+43 1 801 40 0

Belgium
+32 2 482 30 30

France
+33 2 2803 2000

Germany national toll free
08001-536 376

Germany international
+49 6184 90 6940

Italy
+39 02 02 95059 341

Netherlands
+31 76 571 4440

Nordic countries
+358 9 329 100

Russia / CIS
+7 (812) 703 42 15

Spain / Portugal
+34 93 223 09 18

Switzerland
+41 44 454 12 12

UK / Ireland
+44 870 609 9203

Asia:
China
+86 21 6865 4588 or
+86 10 8419 3588

India
+91 22 6716 2200

Japan
+81 45 453 9220

Other Asian countries
+852 2885 4613

Countries not listed:
+49 6184 90 6940 or
+33 2 2803 2000

www.thermo.com/centrifuge

© 2008 Thermo Fisher Scientific Inc. All rights reserved. PLUS is a trademark of Invitrogen. Lipofectamine is a trademark of Life Technologies, Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.