

# Certificate of Analysis

## Product Information:

<b>Orion Number</b>	181811
<b>Description</b>	!REAGENT W/TUBING 3XL BTLS
<b>Lot Number</b>	HR1
<b>Expiration Date</b>	2 Years after Receipt of Product
<b>Method of Analysis</b>	Potentiometric
<b>Acceptance Range</b>	13.1 M +/- 0.6 M Acetic Acid

## Actual Analysis:

<u>Analysis</u>	<u>Value</u>	<u>Unit of Measure</u>
1	13.1	Moles/Liter Acetic Acid

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Potassium Hydrogen Phthalate	185	g

**Certificate Date:** 09/25/03

**Approved By:**



**Quality Control Supervisor**