

Certificate of Analysis

Product Information:

Orion Number	181140
Description	SODIUM STANDARDS KIT Standard 1
Lot Number	HP1
Expiration Date	12/05
Method of Analysis	Flame Photometer
Acceptance Range	19.1 ppm +/- 0.5 ppm Sodium Chloride (NaCl) as Sodium (Na+)

Actual Analysis:

<u>Analysis</u>	<u>Results</u>	<u>Unit of Measure</u>
1	19.1	ppm Sodium Chloride as Sodium

Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

Certificate Date: 12/09/03

Approved By:

Edward Henning

Quality Control Supervisor