

Certificate of Analysis

Product Information:

Orion Number 181811
Description !REAGENT W/TUBING 3XL BTLS
Lot Number EY1
Expiration Date 2 Years after Receipt of Product
Method of Analysis Potentiometric

Actual Analysis:

| <u>Analysis</u> | <u>Results</u> | <u>Unit of Measure</u> |
|-----------------|----------------|---------------------------------|
| 1 | 13.1 | Moles/Liter Glacial Acetic Acid |

Traceability:

| <u>N.I.S.T. Chemical</u> | <u>SRM #</u> | <u>Sublot</u> |
|------------------------------|--------------|---------------|
| Potassium Hydrogen Phthalate | 185 | g |

Certificate Date: 02/23/00

Approved By:



Quality Control Supervisor

Environmental Instruments
Water Analysis

166 Cummings Center

Beverly, MA
01915

(978) 232-6015 Dom Fax
(978) 232-6031 Int'l Fax

www.thermo.com
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

245967-001 Rev. B

Thermo
ELECTRON CORPORATION