

# Certificate of Analysis

## Product Information:

**Orion Number** 9191860  
**Description** ORION PH 9.18 BUFFER PACKOF 5 X 60 ML BOTTLES  
**Lot Number** JY1  
**Expiration Date** 2-11-07  
**Method of Analysis** Potentiometric

## Actual Analysis:

<u>Analysis</u>	<u>Specification</u>	<u>Results</u>	<u>Unit of Measure</u>
1	9.176- 9.190	9.182	pH

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Disodium Hydrogen Phosphate	186 II	g
Potassium Dihydrogen Phosphate	186 I	g
Sodium Bicarbonate	191	C
Sodium Carbonate	191	C

**Certificate Date:** 04/12/05

**Approved By:**



**Quality Control Supervisor**

Environmental Instruments  
Water Analysis

166 Cummings Center

Beverly, MA  
01915

(978) 232-6015 Dom Fax  
(978) 232-6031 Int'l Fax

www.thermo.com  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION

245968-001 Rev. B