

# Certificate of Analysis

## Product Information:

**Orion Number** 011006  
**Description** COND/TDS STD 12.9MS/CM 5x60ML  
**Lot Number** JY1  
**Expiration Date** 2/07  
**Method of Analysis** Potentiometric

## Actual Analysis:

<u>Analysis</u>	<u>Specification</u>	<u>Results</u>	<u>Unit of Measure</u>
1		0.1001	Moles/Liter Potassium Chloride
2	12.9mS/cm+/-0.07mS/cm	12.90	mS/cm at 25 deg. C

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA
Potassium Chloride	2202	NA

**Certificate Date:** 02/24/05

**Approved By:**



**Quality Control Supervisor**

Environmental Instruments  
Water Analysis

166 Cummings Center

Beverly, MA  
01915

(978) 232-6015 Dom Fax  
(978) 232-6031 Int'l Fax

www.thermo.com  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

245968-001 Rev. B

**Thermo**  
ELECTRON CORPORATION