

# Certificate of Analysis

## Product Information:

**Orion Number** 150071  
**Description** REF. ELEC FILL SOLN 5 x 60 ML  
**Lot Number** NNV1  
**Expiration Date** N/A; Recommendation usage 2 years from date of receipt  
**Method of Analysis** Potentiometric

## Actual Analysis:


<u>Analysis</u>	<u>Specification</u>	<u>Results</u>	<u>Unit of Measure</u>
1	2.00 +/- 0.05 M/L	1.99	Moles/Liter Potassium Chloride

## Traceability:

<u>N.I.S.T. Chemical</u>	<u>SRM #</u>	<u>Sublot</u>
Sodium Chloride	2201	NA

**Certificate Date:** 05/29/09

**Approved By:**



**Quality Representative**

**Environmental Instruments**  
Water Analysis Instruments

166 Cummings Center

Beverly, MA  
01915

(978) 232-6000  
(978) 232-6015 fax

[www.thermo.com/water](http://www.thermo.com/water)  
(800) 225-1480 Tech Info

Analyze • Detect • Measure • Control™

245968-001 Rev. D

**Thermo**  
ELECTRON CORPORATION