

SEALING & STAKING

Rapid & Reliable



THE PERFECT PARTNER FOR YOUR SEALER

Combine the industry's best-in-class solutions for fast, dependable protection and handling of microplates in your laboratory. The Thermo Scientific RapidStak™ stacker and our ABgene ALPS® automated sealers were made for each other. They liberate lab personnel for important discovery work with potent productivity and easy, walkaway automation capability.



SEAL & STAK

**FAST,
RELIABLE,
LIBERATING**



The RapidStak microplate stacker and ABgene ALPS automated sealers — both from Thermo Fisher Scientific — make perfect partners.

The Thermo Scientific RapidStak is one of the few instrument loaders that matches the blazing speed and reliable operation of our ABgene ALPS sealers. Simply load the RapidStak with plates — and walk away with confidence.

Perform setup in minutes. No tools required!

- **Fast** — seals up to 210 plates per hour
- **Reliable** — processes up to 150-15 mm plates per run with RapidStak 2X model
- **Easy** — RapidStak and ALPS integrate seamlessly
- **User-friendly** — simple setup; no PC or scheduling software needed
- **Dependable** — backed by the Thermo Fisher Scientific reputation for proven automation solutions

Ideal Plate Preparation Solution

- PCR/QPCR
- Sequencing reactions
- Sample storage
- Colorimetric/fluorometric assays

Easy Upgrade Option

The RapidStak is one of the few stackers that can automate any two instruments simultaneously. Just use our Polara RS software to add one of our Multidrop® automated dispensers. Result: an almost endless array of dispense and seal applications.



Part No.	Product
F01668	RapidStak, Microplate Stacker
F01669	RapidStak 2x, Microplate Stacker
F01362	30 Plate Stak, V3
F01363	50 Plate Stak, V3
AB3000	ALPS 3000 Heat Sealing Unit

© 2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

BRO-SEALSTAK-USA-6/07 BRO-LAI-S&S-6/07 0717011

LIFE SCIENCE RESEARCH
Information, Sales, and Support

3747 N. Meridian Road
P.O. Box 117
Rockford, IL 61105
E-mail: sealandstak@thermo.com
P – 800-445-2812
F – 585-654-4810

Thermo
SCIENTIFIC