














Matrix Tips Specifications



Matrix Tip	Part No.	Finished Good No.	Tip Diameter	Tip Shoulder	Overall Length	Inner Orifice Diameter
 P10 Tip and Filter Tip	18720	7610, 7611, 7612, 7615	0.232~/5.89mm	0.156~/3.96mm	1.229~/31.21mm	0.016~/0.41mm
 20µL Tip and Filter Tip	19563	7620, 7621, 7622, 7625	0.266~/6.76mm	0.141~/3.58mm	1.795~/45.59mm	0.017~/0.431mm
 30µL TallTip and Filter TallTip	19564	7600, 7631, 7632, 7635	0.216~/5.49mm	0.265~/6.73mm	1.60~/40.6mm	0.016~/0.41mm
 250µL Tip	19559	7250, 7151, 7152, 7155, 7156	0.295~/7.49mm	0.635~/16.13mm	1.94~/49.2mm	0.018~/0.46mm
 200µL Filter Tip	19743	7275	0.280~/7.11mm	0.732~/18.59mm	2.13~/54.1mm	0.020~/0.51mm
 250µL TallTip and Filter Tip	19226	7280, 7281, 7282, 7285	0.295~/7.49mm	0.750~/19.05mm	3.06~/77.7mm	0.018~/0.46mm
 300µL Tip and Filter Tip	19635	7320, 7321, 7322, 7325	0.275~/6.95mm	0.750~/19.05mm	2.28~/57.8mm	0.018~/0.46mm
 300µL TallTip and Filter Tip	19227	7080, 7081, 7082, 7085	0.295~/7.49mm	0.750~/19.05mm	4.01~/101.9mm	0.018~/0.46mm
 1250 µL Tip and Filter Tip	19636	8050, 8041, 8042, 8046, 8045	0.348~/8.84mm	1.140~/28.95mm	3.10~/78.7mm	0.035~/0.89mm
 1250µL TallTip and Filter Tip	22635	8250, 8241, 8242, 8246, 8245	0.348~/8.84mm	1.140~/28.95mm	4.03~/102.4mm	0.030~/0.76mm
 12.5µL 384 Tip and Filter Tip	19352	7421, 7422, 7425	0.170~/4.32mm	0.354~/8.99mm	1.60~/40.6mm	0.010~/0.25mm
 30µL 384 Tip and Filter Tip	19353	7431, 7432, 7435	0.170~/4.32mm	0.354~/8.99mm	1.90~/48.3mm	0.012~/0.30mm
 125µL 384 Tip and Filter Tip	19354	7441, 7442, 7445	0.170~/4.32mm	0.354~/8.99mm	2.14~/54.3mm	0.018~/0.46mm

Note: photos not to exact scale

US: Phone: 800 345 0206
 Europe: Phone: +44 (0) 161 486 2110
 Asia: Phone: +81 (0) 47 355 6491

e-mail: info@matrixtechcorp.com
 e-mail: info@matrixtechcorp.eu.com
 e-mail: asia@matrixtechcorp.com

MATRIX
 fluid thinking.
www.matrixtechcorp.com