

AH-629 (17 ml) Swinging Bucket Aluminum Rotor

Rotor Specifications	
Capacity	6 x 17
Tube Dimensions \varnothing x L (mm)	16 x 102
Angle $^{\circ}$	90
Maximum Speed (rpm)	29,000
K Factor	284
Net Weight (kg)	6.2

RCF/Radius	RCF (g)	Radius (cm)
Maximum	155,846	16.59
Average	108,219	11.52
Minimum	60,591	6.45



AH-629 for 17 ml Tubes Rotor Package

Cat. No.	Description
54285	AH-629 for 17 ml Tubes Rotor Package Includes:
51942	Rotor Stand
65937	Vacuum Grease
61556	Lubricant
52030	17 ml Buckets (6)
51362	Extra Overspeed Decal, 29,000 rpm
63712	Bucket O-ring (6)
52240	Bucket Rack
52384	Ultracentrifuge Rotor Log Book
54276	Instruction Manual

Accessories

Cat. No.	Description
52035	36 ml Buckets (6)
03920	Ultracrimp/Clearcrimp Sealing Tool
03919	Crimp Gauge Replacement
03046	Extra Plugs and Caps
12826	Tube Removal Tool
03327	Tube Racks (2)
4209	Conical Adapter Removal Tool (Parts only available directly from Seton Scientific)

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. \varnothing x L (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:
							Cat. No.	Qty/ pkg.	Description	Cat. No.	Qty/ pkg.	No. of places	
06750	17	17.9	PET Thin-Walled Tube	50	29,000	16 x 102	-	-	-	-	-	-	-
03126	17	16	PA Thin-Walled Tube	50	29,000	16 x 102	-	-	-	-	-	-	-
03409	15.5	15.5	PA Ultracrimp Tube	50	29,000	16 x 100	Includes	25	Plug and Cap	-	-	-	03327 Tube Racks, 03920 Crimp Sealing Tool
7066	13	12.7	PET Ultracone Tube (Parts only available directly from Seton Scientific. For use <10C only)	N/A	29,000	16 x 96				4229	1	1	4209 Conical Adapter Removal Tool
5066	13	12.3	PA Ultracone Tube (Parts only available directly from Seton Scientific)	N/A	29,000	16 x 96				4229	1	1	4209 Conical Adapter Removal Tool