

A-621 Fixed Angle Aluminium Rotor

Rotor Specifications	
Capacity (ml)	6 x 250
Tube Dimensions \varnothing x L (mm)	61 x 153
Angle $^{\circ}$	23
Maximum Speed (rpm)	21,000
K Factor	751
Net Weight (kg)	13.64
RCF/Radius	RCF (g)
Maximum	66,993
Average	42,536
Minimum	18,078
RCF/Radius	Radius (cm)
Maximum	13.60
Average	8.635
Minimum	3.67



A-621 Rotor Package

Cat. No.	Description
11876	A-621 Rotor Package Includes:
11998	Rotor Cover
11999	Locking Nut
51942	Rotor Stand
65937	Vacuum Grease
61556	Lubricant
52277	Extra Overspeed Decal, 21,000 rpm
92233	Rotor Cover O-ring, 6 Extra
90793	Locking Nut O-ring, 6 Extra
52384	Ultracentrifuge Rotor Log Book
11872	Instruction Manual
52392	Rotor Cover Tool
65614	Retaining Ring Pliers

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. \varnothing x L (mm)	Sealing Assembly Required:		Accessories Required:	
							Cat. No.	Qty/ pkg		
314250 H01	250	200	PA Oak Ridge Bottle	4	21,000	61 x 120	314249 H01	1	Delrin Screw Cap	-
03069	250	280	PP Dry-Spin Bottles	12	13,000	61 x 146	Included	12	Sealed Lid	-
03098	250	270	PC Dry-Spin Bottles	12	16,000	61 x 146	Included	12	Sealed Lid	-