

Stepsaver 65V13 Vertical Titanium Rotor

Rotor Specifications		
Capacity (ml)	8 x 6	8 x 13.5
Tube Dimensions \varnothing x L (mm)	13 x 66	16 x 93
Angle $^{\circ}$	0	0
Maximum Speed (rpm)	65,000	65,000
K Factor	10.3	12.9
Net Weight (kg)	9.7	9.7

RCF/Radius	RCF (g)	R (cm)	RCF (g)	R (cm)
Maximum	393,120	8.33	400,671	8.49
Average	361,972	7.67	361,972	7.67
Minimum	330,825	7.01	323,274	6.85



Stepsaver 65V13 Rotor Packages

Cat. No.	Description
79243	Stepsaver 65V13 for 13.5 ml Tubes Rotor Package Includes:
51942	Rotor Stand
51345	Extra Overspeed Decal, 65,000 rpm
52384	Ultracentrifuge Rotor Log Book
79231	Instruction Manual
79246	PA Ultracrimp Tube
79296	13.5 ml Ultracrimp Cartridge Set (8)
79247	Tube Racks (2)
79212	Cartridge Set Rack
79285	Stepsaver 65V13 for 13.5 ml Beckman Quick-Seal Tubes Complete Rotor Package Includes:
51942	Rotor Stand
51345	Extra Overspeed Decal, 65,000 rpm
52384	Ultracentrifuge Rotor Log Book
79231	Instruction Manual
79298	13.5 ml Quick-Seal Cartridge Set (8)
79212	Cartridge Set Rack
79290	Stepsaver 65V13 for 6 ml Tubes Rotor Package Includes:
51942	Rotor Stand
51345	Extra Overspeed Decal, 65,000 rpm
52384	Ultracentrifuge Rotor Log Book
79231	Instruction Manual
79273	PA Ultracrimp Tube
79297	6 ml Ultracrimp Cartridge Set (8)
03921	Tube Racks (2)
79212	Cartridge Set Rack

Accessories

Cat. No.	Description
03920	Ultracrimp/Clearcrimp Sealing Tool
03919	Crimp Gauge Replacement
79296	13.5 ml Ultracrimp Cartridge Set (8)
79244	13.5 ml Ultracrimp Individual Cartridge (1)
79297	6 ml Ultracrimp Cartridge Set (8)
79260	6 ml Ultracrimp Individual Cartridge (1)
79298	13.5 ml Quick-Seal Cartridge Set (8)
79252	13.5 ml Quick-Seal Individual Cartridge (1)
79212	Cartridge Set Rack

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. Ø x L (mm)	Sealing Assembly Required:			Accessories Required:
							Cat. No.	Qty/ pkg	Description	
79246	13.5	13.8	PA Ultracrimp Tube	50	65,000	16 x 93	Included	50	Plug and Cap	79247 Tube Racks, 03920 Crimp Sealing Tools, 79296 Cartridges
79294	13.5	13.8	PET Clearcrimp Tube	50	65,000	16 x 93	Included	50	Plug and Cap	79247 Tube Racks, 03920 Crimp Sealing Tools, 79296 Cartridges
342413	13.5	13.5	PA Quick-Seal Tube (Tubes only available from Beckman)	-	65,000	16 x 76	-	-	-	79298 Cartridges
344322	13.5	13.5	PET Quick-Seal Tube (Tubes only available from Beckman)	-	65,000	16 x 76	-	-	-	79298 Cartridges
79273	6	5.7	PA Ultracrimp Tube	50	65,000	13 x 66	Included	50	Plug and Cap	03921 Tube Racks, 03920 Crimp Sealing Tool, 79297 Cartridges
79293	6	5.7	PET Clearcrimp Tube	50	65,000	13 x 66	Included	50	Plug and Cap	03921 Tube Racks, 03920 Crimp Sealing Tool, 79297 Cartridges