

Surespin 630 (17 ml) Swinging Bucket Titanium Rotor

Rotor Specifications	
Capacity (ml)	6 x 17
Tube Dimensions \varnothing x L mm	16 x 102
Angle $^{\circ}$	90
Maximum Speed (rpm)	30,000
K Factor	268
Net Weight (kg)	8
RCF/Radius	RCF (g)
Maximum	166,880
Average	115,610
Minimum	64,339
RCF/Radius	Radius (cm)
Maximum	16.60
Average	11.5
Minimum	6.4



Surespin 630 for 17 ml Tubes Rotor Package

Cat. No.	Description
79363	Surespin 630 (17 ml) Rotor Package Includes:
51942	Rotor Stand
65937	Vacuum Grease
61556	Lubricant
79364	17 ml Buckets (6)
51361	Extra Overspeed Decal, 30,000 rpm
79370	Bucket O-ring (6)
52240	Bucket Rack
79369	Instruction Manual
52384	Ultracentrifuge Rotor Log Book

Accessories

Cat. No.	Description
79365	36 ml Buckets (6)
03920	Ultracrimp/Clearcrimp Sealing Tool
03919	Crimp Gauge Replacement
03046	Extra Plugs and Caps
12826	Tube Removal Tool
03327	Tube Racks (2)
4209	Conical Adapter Removal Tool (Parts only available directly from Seton Scientific)

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. \varnothing x L (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:
							Cat. No.	Qty/ pkg.	Description	Cat. No.	Qty/ pkg.	No. of places	
06750	17	17.9	PET Thin-Walled Tube	50	30,000	16 x 102	-	-	-	-	-	-	-
03126	17	16	PA Thin-Walled Tube	50	30,000	16 x 102	-	-	-	-	-	-	-
03409	15.5	15.5	PA Ultracrimp Tube	50	30,000	16 x 100	Included	50	Plug and Cap	-	-	-	03327 Tube Racks, 03920 Crimp Sealing Tool
7066	13	12.7	PET Ultracone Tube (Parts only available directly from Seton Scientific. For use <10C only)	N/A	30,000	16 x 96	-	-	-	4229	1	1	4209 Conical Adapter Removal Tool
5066	13	12.3	PA Ultracone Tube (Parts only available directly from Seton Scientific)	N/A	30,000	16 x 96	-	-	-	4229	1	1	4209 Conical Adapter Removal Tool