

TFT-80.2 Fixed Angle Titanium Rotor

Rotor Specifications	
Capacity (ml)	12 x 2
Tube Dimensions \varnothing x L (mm)	8 x 51
Angle $^{\circ}$	25
Maximum Speed (rpm)	80,000
K Factor	22.6
Net Weight (kg)	2.5
RCF/Radius	RCF (g)
Maximum	429,643
Average	335,994
Minimum	242,344
RCF/Radius	Radius (cm)
Maximum	6.01
Average	4.7
Minimum	3.39



TFT-80.2 Rotor Package

Cat. No.	Description
54456	TFT-80.2 Rotor Package Includes:
54463	Rotor Cover
54464	Locking Nut
51942	Rotor Stand
65937	Vacuum Grease
61556	Lubricant
11559	Rotor Cover Tool
11628	Extra Overspeed Decal, 80,000 rpm
54465	Rotor Cover O-ring, 5 Extra
54466	Locking Nut O-ring, 5 Extra
52384	Ultracentrifuge Rotor Log Book
54462	Instruction Manual

Accessories

Cat. No.	Description
11538	Tube Removal Tool

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. \varnothing x L (mm)	Sealing Assembly Required:		Accessories Required:	
							Cat. No.	Qty/ pkg		
54460	2	1.3	PA Thin-Walled Tube	50	80,000	8 x 51	54459	1	Multipiece Aluminum Assembly	11538 Tube Removal Tool