

TH-660 Swinging Bucket Titanium Rotor

Rotor Specifications	
Capacity (ml)	6 x 4.4
Tube Dimensions \varnothing x L (mm)	11 x 60
Angle $^{\circ}$	90
Maximum Speed (rpm)	60,000
K Factor	44.4
Net Weight (kg)	5.8

RCF/Radius	RCF (g)	Radius (cm)
Maximum	488,576	12.15
Average	374,173	9.305
Minimum	259,770	6.46



TH-660 Rotor Package

Cat. No.	Description
10502	TH-660 Rotor Package Includes:
51942	Rotor Stand
65937	Vacuum Grease
61556	Lubricant
10509	4.4 ml Bucket Set (6)
51346	Extra Overspeed Decal, 60,000 rpm
63713	Bucket O-ring (6)
52243	Bucket Rack
73031	Bucket Cap Tool
10506	Instruction Manual

Accessories

Cat. No.	Description
4209	Conical Adapter Removal Tool (Parts only available directly from Seton Scientific)
45537	Re-Seal 11 mm x 32 mm Tube Vise

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. \varnothing x L (mm)	Sealing Assembly Required:			Adapters Required:			Accessories Required:
							Cat. No.	Qty/ pkg	Description	Cat. No.	Qty/ pkg	No. of places	
03955	4.4	4	PA Thin-Walled Tube	50	60,000	11 x 60	-	-	-	-	-	-	-
7062	3.4	3.4	PET Thin-Walled Ultracone Tube (Parts only available directly from Seton Scientific. For use <10C only)	N/A	60,000	11 x 57	-	-	-	4226	1	-	4209 Conical Adapter Removal Tool
5062	3.4	3.4	PA Thin-Walled Ultracone Tube (Parts only available directly from Seton Scientific)	N/A	60,000	11 x 57	-	-	-	4226	1	-	4209 Conical Adapter Removal Tool
9012	2	2	PA Cone-Top Tube (Parts only available directly from Seton Scientific)	N/A	60,000	11 x 32	4246	1	11mm Crown	4075	10	-	4088 Wrench, 4250 Tube Removal Tool, 4172 Tube Vise