

TV-860 Vertical Titanium Rotor

Rotor Specifications	
Capacity (ml)	8 x 35
Tube Dimensions \varnothing x L (mm)	25 x 97
Angle $^{\circ}$	0
Maximum Speed (rpm)	60,000
K Factor	25.2
Net Weight (kg)	11.4
RCF/Radius	RCF (g)
Maximum	340,596
Average	289,325
Minimum	238,055
RCF/Radius	Radius (cm)
Maximum	8.47
Average	7.195
Minimum	5.92



TV-860 Rotor Package

Cat. No.	Description
53289	TV-860 Rotor Package Includes:
51942	Rotor Stand
61556	Lubricant
51346	Extra Overspeed Decal, 60,000 rpm
52384	Ultracentrifuge Rotor Log Book
53292	Instruction Manual

Accessories

Cat. No.	Description
03920	Ultracrimp/Clearcrimp Sealing Tool
03919	Crimp Gauge Replacement
03999	Extra Plugs and Caps
12771	Tool Kit for Ultracrimp/Clearcrimp Tubes (Includes: Rotor Holding Fixture (52540), Torque Wrench (52580), Socket (54519))

Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty	Max. Speed (rpm)	Dims. \varnothing x L (mm)	Sealing Assembly Required:			Accessories Required:
							Cat. No.	Qty/ pkg	Description	
03989	35	35	PA Ultracrimp Tube	25	60,000	25 x 97	Included	50	Plug and Cap	03994 Rotor Cap, 03935 Tube Racks, 03920 Crimp Sealing Tool, 12771 Tool Kit