

# TV-1665 Vertical Titanium Rotor

Rotor Specifications	
Capacity (ml)	16 x 6
Tube Dimensions $\varnothing$ x L (mm)	13 x 66
Angle $^{\circ}$	0
Maximum Speed (rpm)	65,000
K Factor	9.8
Net Weight (kg)	9.07

  

RCF/Radius	RCF (g)	Radius (cm)
Maximum	414,357	8.78
Average	383,209	8.12
Minimum	352,062	7.46



## TV-1665 Rotor Package

Cat. No.	Description
<b>11672</b>	<b>TV-1665 Rotor Package Includes:</b>
51942	Rotor Stand
61556	Lubricant
51345	Extra Overspeed Decal, 65,000 rpm
52384	Ultracentrifuge Rotor Log Book
11670	Instruction Manual

## Accessories

Cat. No.	Description
03920	Ultracrimp/Clearcrimp Sealing Tool
03919	Crimp Gauge Replacement
03999	ULTRACRIMP Plugs and Caps (25)
12770	Tool Kit for Ultracrimp Sealing System (Includes: Rotor Holding Fixture (52540), Torque Wrench (52580), Socket Cap Wrench (52882))
52598	Tool Kit for 5 ml Thin-Walled Tube (03127) (Includes: Wrench (52363), Rotor Holding Fixture (52540))

## Tube Ordering Information

Cat. No.	Tube Vol. (ml)	Vol. Fill (ml)	Description	Qty.	Max. Speed (rpm)	Dims. $\varnothing$ x L (mm)	Sealing Assembly Required:			Accessories Required:
							Cat. No.	Qty/ pkg	Description	
03945	6	5.7	PA Ultracrimp Tube	50	65,000	13 x 66	Included	50	Plug and Cap	03991 Rotor Cap, 03921 Tube Racks, 03920 Crimp Sealing Tool, 12770 Tool Kit
01086	6	5.7	PET Clearcrimp Tube	50	65,000	13 x 66	Included	50	Plug and Cap	03991 Rotor Cap, 03921 Tube Racks, 03920 Crimp Sealing Tool, 12770 Tool Kit
03127	5	4.4	PA Thin Walled Tubes	50	65,000	13 x 51		1	Body Cap	52598 Tool Kit
							52521	8	Plugs (Body Cap and Plugs must be used together)	