

TZ-28 Zonal Titanium Rotor

Rotor Specifications	
Capacity (ml)	1350
Tube Dimensions \varnothing x L (mm)	N/A
Angle $^{\circ}$	90
Maximum Speed (rpm)	28,000
K Factor	312
Net Weight (kg)	7.8

RCF/Radius	RCF (g)	Radius (cm)
Maximum	83,547	9.53
Average	57,579	6.575
Minimum	31,701	3.62



TZ-28 Zonal Titanium Rotor Packages:

Cat. No.	Description
52358	TZ-28 Zonal Titanium Rotor Includes:
49609	Core Septa
66193	O-ring for Septa (12)
49724	Rotor Cover Wrench
49026	T-wrench
52540	Rotor Holding Fixture
66026	Sample Inject Needle
61556	Lubricant
65755	O-ring for Rotor Cover (2)
49007	Feed Adapter
63037	O-ring for Feed Adapter
12271	TZ-28 Continuous Flow Ultracentrifuge Conversion Kit Includes (required):
65937	Vacuum Grease
52384	Ultracentrifuge Rotor Log Book
52415	Instruction Manual
52496	Tapered Drive Assembly
52418	Rotor Adapter for Drive Assembly
60912	Socket Head Screw for Rotor Adapters (3)
52497	Drive Adapter Assembly Wrench
51363	Overspeed Decal, 28,000 rpm
52417	Sealing Cover
66330	O-ring for Sealing cover (2)
62987	Socket Head Screw for Sealing Cover (3)
52425	Distributor Body
64578	O-ring for Upper Distributor Body (2)
66308	O-ring for Lower Distributor Body (2)
52416	Hold Down Plate for Distributor Body
61289	O-ring for Hold Down Plate (12)
66329	Socket Head Screw for Hold Down Plate (3)