

Thermo Scientific ABsolute™ QPCR Seal

Description: The ABsolute™ QPCR Seal is an optically clear, non-tacky pressure-sensitive adhesive developed specifically for QPCR and other fluorescence applications.

Ordering Information: AB-1170 ABsolute™ QPCR Seal 50 sheets

Sheet size: Length: 141.3 ± 1.6mm (with tabs)
120.7 ± 1.6mm (without tabs)
Width: 79.4 ± 1.6mm

Storage: Store in a cool, dry place. Avoid direct exposure to sunlight. It is recommended to use these seals within one year from date of purchase.

Physical Properties: Seal thickness: 100µm (without release liner)
Temperature range: -80°C to 110°C

Sealing: AB-1170: Bond strength is dependent on application pressure, so sufficient pressure must be applied. For optimal seal application, ABgene® recommend the use of the Adhesive Seal Applicator (AB-1391), available free of charge when ordered with AB-1170.

Troubleshooting

For technical information or troubleshooting contact Thermo Scientific Genomics Tech Support:

| | United States | Rest of World |
|--------------|--|--|
| Email | techservice.genomics@thermofisher.com | techservice.emea.genomics@thermofisher.com |
| Phone | +1 (800) 235-9880 | +44 (0)1372 840410 |

ABsolute™ and SEAL-IT 100™ are trademarks of Thermo Fisher Scientific and its subsidiaries

© 2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Literature Code: AB-1170

For Research Purposes Only