

Adhesive Crystallography Seal

Description: An optically clear, non-tacky adhesive plate seal ideal for crystallography, PCR and genotyping applications.

Ordering Information: AB-1450 Adhesive Crystallography Seal Adhesive 50 sheets
AB-2450 Crystallography Seal Roll for SEAL-IT 250m roll
100™

Sheet size: Length: 141.3 ± 1.6mm (with tabs)
(AB-1450 only) 120.7 ± 1.6mm (without tabs)
Width: 79.4 ± 1.6mm

Storage: Store in a cool, dry place. Avoid direct exposure to sunlight. It is recommended to use these seals within one year from date of purchase.

Physical Properties: Seal thickness: 100µm (without release liner)
Temperature range: -70°C to 110°C

Sealing: AB-1450: Bond strength is dependent on application pressure, so sufficient pressure must be applied. For optimal seal application, Thermo Fisher Scientific recommend the use of the Adhesive Seal Applicator (AB-1391), available free of charge when ordered with AB-1450.

AB-2450: For recommended sealing conditions please refer to the SEAL-IT 100™ Instruction Manual.

SEAL-IT 100™ is a trademark of Thermo Fisher Scientific and its subsidiaries

For Research Purposes Only