

AutoSCAN™ Software Training Course Outline

	Module 1	Module 2	Module 3	Lunch
Session	8:00 - 8.30 am	8:30 – 9:30 am	9:30 - 11:30 am	11:30-12:30 pm
Day 1	Course Basics Welcome Introductions Get Acquainted Course Timeline Materials Overview Course Objectives Detailed Course Outline	Product Information Product Definition Compatible Thermo Products SCADA Components Program Capabilities: What is Poller What is Viewer What is Vox What is EXPORT Typical AutoSCAN Setup Communication Options System Requirements/ Installation Procedures Desktop Shortcuts Test Your Knowledge Exercise	Poller (Hands-On) Review of Poller's Capabilities Starting /Initializing Poller Steps to Setup: System Parameter Setup Driver Parameter Setup Unit Scan Parameter Setup Manual Scan Cancelling a Scan Debugging System Error Log Log Drivers Creating a Backdoor Protocol Manual Import Unit Initialization	Approximate Time
	Module 4	See next page for day 2		
Session	12:30 - 5:00 pm			
Day 1 (Cont)	Viewer (Hands-On) Starting Viewer Initializing Viewer: Test Your Knowledge (Initializing Viewer) Search RTU Icon: Test Your Knowledge (Search RTU) Instantaneous Data Icon: Test Your Knowledge (Instantaneous Data Icon) Static Data Icon: Test Your Knowledge (Changing Static Data)			

	Module 4 (Cont)	Lunch	Module 4 (Cont)	Module 5
Session	8:00 – 11:30 am	11:30 – 12:30 pm	12:30 - 2:30 pm	2:30-5:00 pm
Day 2	<p>Viewer (Hands-On) System Data Icon: Test Your Knowledge (System Data) Historical Daily Data Icon: Test Your Knowledge (Historical Daily Data) Historical Periodic Data Icon: Test Your Knowledge (Historical Periodic Data) Historical Alarm Data Icon: Test Your Knowledge (Historical Alarm Data) Historical Audit Trail Data Icon: Test Your Knowledge (Historical Audit Data) Print History Icon: Test Your Knowledge (Print History Icon) Summary Report Icon: Test Your Knowledge (Summary Report Icon) Alarm Summary Icon: Test Your Knowledge (Alarm Summary Report) Unacknowledged Alarms</p>	<p>Approximate Time</p>	<p>Viewer (Hands-On) Calculator Icon Test Your Knowledge (Calculator) Enable/Disable Function Icon: Inputs Outputs Run Switching Gas Sampling PID Control Chromatograph History Function: Display Edit Tank Gauging Test Your Knowledge (Enable/Disable Icon) Manual Scan Icon Test Your Knowledge (Manual Scan Icon) Unacknowledged Alarms Icon Graph Historical Data Icon Test Your Knowledge (Graph Historical Data) Reports: Weekly/Daily Management Daily Management Report Configurations: Summary Report Type: Basic Basic w/Tubing & Casing Daily Report Type: Basic Water Flow-Cal Type Location Station</p>	<p>VOX Starting VOX Main Menu File Functions: Printer Setup Data Directory Exit Configuration Functions: Configure Call-In Data Window: Cascade Tile: Horizontal Vertical Always on Top Toolbar: Data Directory Configure Call-Out Data View/Record Voice File Configure Call-IN Data Create Database Phone Book Call-In Status Call-Out Status Log Setup Local Alarm Alert Type Alarm Log Convert VOX to WAV Ack Call-Out Alarm Lock/Unlock Print Function Touch Tone Phone Operation: Call-In Operation Call-Out Operation</p>