

NOT TO SCALE

127.76 ±0.25

103.50 ±0.25

9.50 ±0.25

12.13 ±0.25

4.50

4.00 (X8)

∅3.00
±0.10

0.40

9.00 ±0.10

85.48 ±0.25

67.50 ±0.25

8.99 ±0.25

4.50

38.00

38.00

30.00

30.00

17.60

2.00 ±0.25

VIAL FORM 384pos

PLATES STACKED

SECTION ON CENTRE LINE

TITLE

Diamond Ultra 384

DRG. NO.

AB-2150

ISSUE

ORIGINAL

DATE

26/5/09

All dimensions in mm

Copyright Thermo Fisher Scientific
This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.