

Pipette Compatibility Chart for P/N: 990-02, 990-03 and 990-04

MANUFACTURER	MODEL	VOLUME (μL)
THERMO SCIENTIFIC FINNPIPETTE	NOVUS ELECTRONIC, SINGLE AND MULTICHANNEL	100-1200
	FOCUS, SINGLE CHANNEL	100-1000
	F1, SINGLE CHANNEL	100-1000
	F2, SINGLE CHANNEL	100-1000
BIOHIT	PROLINE	100-1000
EPPENDORF	PHYSIO CARE RESEARCH, SINGLE CHANNEL	100-1000
	PHYSIO CARE RESEARCH, MULTICHANNEL	50-1200
	2000 RESEARCH	200-1000
	2100 RESEARCH	100-1000
FISHER BRAND	FINNPIPETTE	2000-1000
	FINNPIPETTE II	100-1000
	FINNPIPETTE II	200-1000
GILSON	P-100 PIPETMAN	100-1000
HAMILTON	SOFTGRIP	100-1000
NICHIRYO	5000V-M TRIPLE VOLUME	200-1000
OXFORD	BENCHMATE	200-1000
RAININ	R-1000 PIPET-PLUS	100-1000
SOCOREX	821	200-1000
	825 ACURA	100-1000
	822 CALIBRA	100-1000
TITERTEK	DIGITAL	200-1000
VISTALAB	OVATION	100-1000
VWR	SIGNATURE	100-1000