

Thermo Scientific Abgene Multisip Plugs complement the innovative line of Thermo Scientific storage tubes. The flexibility of individual tube capping combined with the pierceable cap design allows individual sample access even while the storage tubes remain sealed.

Thermo Scientific Abgene® Multisip™ Plugs



Enhanced Materials. Thermo Scientific Abgene Multisip Plugs are made of an advanced thermoplastic elastomer that offers excellent chemical resistance and low vapor transmission rates, minimizing evaporation and water uptake in DMSO-based samples. The plug's V-shaped design reliably guides liquid handling needles to the same location repeatedly, while this elastomer ensures sealing performance will not degrade, even after hundreds of cap piercings.

Application Flexibility. Designed for storage down to the vapor phase of liquid nitrogen, Abgene Multisip Plugs are available in bulk or pre-installed in a wide array of our storage tubes. In addition, optional pre-racked plugs provide the ability to apply 96 plugs at once using the Thermo Scientific Abgene Multisip Applicator.

Knowledgeable Support. Our team has a strong reputation for supporting both equipment and consumables, and can demonstrate how Abgene Multisip Plugs are the ideal solution for your individual tube sealing needs.



Thermo Scientific Abgene Multisip Plugs

Cap Material	Thermoplastic elastomer
Cap Temperature Range	-180 °C to 120 °C
Cap Compatibility	96-format Thermo Scientific Matrix® and Abgene Storage Tubes
Cap Width	7 or 10 mm

Ordering Guide

Item No.	Description	Qty
AB-0972	Abgene Multisip Split Septum Plugs, 7 mm, Bulk	1000 plugs/case
AB-1174	Abgene 96-Format Multisip Septum Plugs, 7 mm	10 trays of 96 plugs/case
AB-0972	Abgene Multisip Split Septum Plugs, 10 mm, Bulk	1000 plugs/case
AB-1227	Abgene Multisip Cap-It Application Tool	Each
AB-1173/1	Abgene Multisip Applicator, includes Abgene Multisip Applicator Adapter for use with Abgene 1.2 ml Tubes	Each
AB-1173/2	Abgene Multisip Applicator, includes Abgene Multisip Applicator Adapter for use with Abgene 0.65 ml Tubes	Each
AB-1302	Abgene Multisip Applicator Adapter for use with Abgene 1.2 ml Tubes	Each
AB-1303	Abgene Multisip Applicator Adapter for use with Abgene 0.65 ml Tubes	Each

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

North America: Tel: 800.345.0206 | email: matrix.info@thermofisher.com
 Europe: Tel: +44 (0) 161 486 2110 | email: matrix.eu.info@thermofisher.com
 Asia: email: matrix.ap.info@thermofisher.com

www.thermo.com/matrix

Thermo
SCIENTIFIC