

Anchored Oligo dT Primers

Description: Thermo Scientific Anchored Oligo-dT Primers are specifically designed for synthesising long cDNAs from mRNAs carrying poly-A tails. Since this primer only anneals to the poly-A tail junction of the mRNA template, it eliminates transcription through the poly-A tail providing more effective cDNA synthesis.

Ordering Information: AB-1247 Anchored Oligo dT Primers (500ng/µl) 100µl

Storage: Store at -20°C until ready for use for up to 1 year. Shipped on ice within the UK and on dry ice for international and within the US.

For Research Purposes Only