

Value Sheet / Datenblatt / Fiche de valeurs / Tabulka hodnot / Φύλλο τιμών / Hoja de valores / Väärtuste leht / Értéklap / Scheda dei valori / Verçiu lapas / Vērtību tabula / Arkusz wartości / Folha de valores / Fișă de valori / Hárok hodnôt / Referensblad

## Myoglobin Calibration Set

REF 981849 4 x 1 ml

LOT G631 2012-11

THIS VALUE SHEET IS APPLICABLE FOR KONELAB™ AND T SERIES OUTSIDE THE US.<sup>1)</sup>

REF	TEST	CAL	⊕	UNIT
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 1	65	μg/l ng/ml
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 2	130	μg/l ng/ml
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 3	260	μg/l ng/ml
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 4	520	μg/l ng/ml

Notes / Hinweis / Remarque / Poznámky / Σημειώσεις / Notas / Märkused / Megjegyzések / Note / Pastabos / Piezīmes / Uwaga / Notas / Note / Poznámky / Kommentarer:

1) This value sheet is applicable for Konelab™ and T Series outside the US. / Das vorliegende Datenblatt gilt für Konelab™ and die T Series außerhalb der USA. / Cette feuille de valeurs s'utilise avec le Konelab™ et la gamme T hors des Etats-Unis. / Tento protokol platí pro přístroje Konelab™ a řady T mimo USA. / Αυτό το φύλλο τιμών εφαρμόζεται για τις σειρές Konelab™ και T Series εκτός των ΗΠΑ. / Esta hoja de valores es aplicable para Konelab™ y T Series fuera de los EE. UU. / Väärtuste tabel on rakendatav Konelab™ ja T Series aparatuuridele väljaspool Ameerika Ühendriike. / Ez az értéktáblázat az Egyesült Államokon kívül forgalmazott Konelab™ és T Seriesre érvényes. / Questo foglio valori è applicabile per Konelab™ e la serie T all'esterno degli USA. / Šis reikšmių lapas taikomas Konelab™ ir „T Series“ už JAV ribų. / Šī vērtību lapa attiecas tikai uz Konelab™ un T Series ārpus ASV. / Niniejsze zestawienie wartości dotyczy Konelab™ oraz T Series poza granicami USA. / Esta folha de valores aplica-se ao Konelab™ e à série T fora dos EUA. / Această fișă de valori se aplică pentru Konelab și Seria T în afara S.U.A. / Tento hárok hodnôt je použiteľný pre Konelab™ a T Sériu mimo USA. / Detta faktablad gäller för Konelab™ och T-serien utanför USA.



Thermo Fisher Scientific Oy  
Clinical Diagnostics Finland  
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland  
Tel. +358 9 329 100, Fax +358 9 3291 0300  
www.thermo.com/konelab  
2011-08-23



Value Sheet / Datenblatt / Fiche de valeurs / Tabulka hodnot / Φύλλο τιμών / Hoja de valores / Väärtuste leht / Értéklap / Scheda dei valori / Verčij lapas / Vērtību tabula / Arkusz wartości / Folha de valores / Fișa de valori / Hárok hodnôt / Referensblad

## Myoglobin Calibration Set

REF 981849 4 x 1 ml

LOT F597 2011-11

THIS VALUE SHEET IS APPLICABLE FOR KONELAB™ AND T SERIES OUTSIDE THE US.<sup>1)</sup>

REF	TEST	CAL	⊕	UNIT
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 1	65	μg/l ng/ml
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 2	130	μg/l ng/ml
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 3	260	μg/l ng/ml
Myoglobin Plus 981923, 981939	Myoglobin Plus	Calibrator 4	520	μg/l ng/ml

Notes / Hinweis / Remarque / Poznámky / Σημειώσεις / Notas / Märkused / Megjegyzések / Note / Pastabos / Piezīmes / Uwaga / Notas / Note / Poznámky / Kommentarer:

1) This value sheet is applicable for Konelab™ and T Series outside the US. / Das vorliegende Datenblatt gilt für Konelab™ and die T Series außerhalb der USA. / Cette feuille de valeurs s'utilise avec le Konelab™ et la gamme T hors des Etats-Unis. / Tento protokol platí pro přístroje Konelab™ a řady T mimo USA. / Αυτό το φύλλο τιμών εφαρμόζεται για τις σειρές Konelab™ και T Series εκτός των ΗΠΑ. / Esta hoja de valores es aplicable para Konelab™ y T Series fuera de los EE. UU. / Väärtuste tabel on rakendatav Konelab™ ja T Series aparatuuridele väljaspool Ameerika Ühendriike. / Ez az értéktáblázat az Egyesült Államokon kívül forgalmazott Konelab™ és T Seriesre érvényes. / Questo foglio valori è applicabile per Konelab™ e la serie T all'esterno degli USA. / Šis reikšmių lapas taikomas Konelab™ ir „T Series“ už JAV ribų. / Šī vērtību lapa attiecas tikai uz Konelab™ un T Series ārpus ASV. / Niniejsze zestawienie wartości dotyczy Konelab™ oraz T Series poza granicami USA. / Esta folha de valores aplica-se ao Konelab™ e à série T fora dos EUA. / Această fișă de valori se aplică pentru Konelab și Seria T în afara S.U.A. / Tento hárok hodnôt je použiteľný pre Konelab™ a T Sériu mimo USA. / Detta faktablad gäller för Konelab™ och T-serien utanför USA.



Thermo Fisher Scientific Oy  
Clinical Diagnostics Finland  
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland  
Tel. +358 9 329 100, Fax +358 9 3291 0300  
www.thermo.com/konelab  
2010-11-18

