

# Material Safety Data Sheet

Not classified as Hazardous according to criteria of Worksafe Australia  
Not classified as Hazardous according to Occupational Safety and Health Administration's (OSHA)  
Communication Standard, 29 CFR 1910.1200

## ALT (GPT)

### 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

THERMO ELECTRON  
189 - 199 Browns Rd  
NOBLE PARK VIC 3174  
AUSTRALIA  
Tel: +61 3 9790 4100  
Fax: +61 3 9790 4155  
E-mail: [info.clinicalchemistry@thermo.com](mailto:info.clinicalchemistry@thermo.com)

THERMO ELECTRON  
331 South 104<sup>th</sup> Street  
LOUISVILLE, CO 80027  
U.S.A  
Tel: (303) 581 6428  
Fax: (303) 581 6429  
E-mail: [info.clinicalchemistry@thermo.com](mailto:info.clinicalchemistry@thermo.com)

**Product Name:** ALT (GPT)

**Catalog Numbers:** 7200-004

**Use:** For in vitro diagnostic use and for the quantitative determination of ALT (GPT) in serum on the Olympus Demand<sup>®</sup> and Reply<sup>®</sup>.

**UN Number:** None allocated

**Proper Shipping Name:** NONE ALLOCATED

**Dangerous Goods Class:** None allocated

**Subsidiary risk:** None allocated

**Packing Group:** None allocated

**Hazchem Code:** None allocated

**Poison Schedule:** None allocated

#### Contact Point

##### Australia

Quality Assurance Manager:

Tel: +61 3 9790 4100

Mon – Fri 9:00am to 5:00pm

##### U.S.A

Chemtel

24 Hour Emergency Assistance

1-800-255-3924

### 2. COMPOSITION / INFORMATION ON INGREDIENTS

SUBSTANCE NAME	Proportion	CAS Number
----------------	------------	------------

All ingredients in this product are not hazardous according to the criteria of Worksafe Australia and OSHA Communication Standard 29 CFR 1910.1200

### 3. HAZARD IDENTIFICATION

Not Classified as hazardous according to the criteria of Worksafe Australia and OSHA Communication Standard 29 CFR 1910.1200

**Hazard Category:** None allocated

#### ACUTE HEALTH EFFECTS

##### Swallowed:

Swallowing large quantities may cause irritation to mouth, throat and stomach, which may lead to nausea, vomiting and diarrhoea.

##### Eye:

May cause mechanical irritation to the eyes, with effects including: tearing, pain, stinging and blurred vision. These effects are anticipated to be of a transient (short acting) nature and no long term injury is envisaged.

##### Skin:

May cause irritation to the skin.

##### Inhaled:

May cause irritation to the mouth, throat and upper respiratory system.

##### Chronic:

Prolonged or repeated inhalation of the dust may cause sneezing, coughing and mucous buildup.

# Material Safety Data Sheet

Not classified as Hazardous according to criteria of Worksafe Australia  
Not classified as Hazardous according to Occupational Safety and Health Administration's (OSHA)  
Communication Standard, 29 CFR 1910.1200

## ALT (GPT)

### 4. FIRST AID MEASURES

**Swallowed:**

If swallowed, **DO NOT** induce vomiting. If conscious, give 1 to 2 glasses of water to drink. Seek immediate medical assistance.

**Eye:**

If material is splashed into eyes, immediately, flush with plenty of water for 15 minutes, ensuring eye lids are held open. If irritation persists transport to hospital or doctor.

**Skin:**

If material is splashed onto the skin, remove any contaminated clothing and wash skin thoroughly with soap and water. If irritation persists transport to hospital or doctor.

**Inhaled:**

Move victim to fresh air. Apply resuscitation if victim is not breathing.

**First Aid Facilities:**

Eye wash fountain, safety shower and normal wash room facilities.

**Advice to Doctor:**

Treat symptomatically.

In case of poisoning, contact Poisons Information Centre

**In Australia call Tel: 131126**

**In New Zealand Tel: 034747000**

### 5. FIRE-FIGHTING MEASURES

**Fire/Explosion Hazard**

If safe to do so, move undamaged containers from fire area.

**Hazardous Decomposition Products:** Decomposes on heating emitting oxides of carbon and noxious smoke.

**Fire Fighting Procedures:** Fire fighters to wear Self-contained breathing apparatus (SCBA) in confined spaces, in oxygen deficient atmospheres or if exposed to products of decomposition. Full protective clothing is also recommended.

**Extinguishing Media:** Use extinguishing media suitable for surrounding fire situation.

**Flammability**

This material is not a combustible or flammable solid.

### 6. ACCIDENTAL RELEASE MEASURES

Avoid generating dusts. Wear suitable protective equipment. Ventilate area. If possible wet area down to prevent high dust levels. If available, use dustless methods, such as a HEPA vacuum and filter. Otherwise, use a non-sparking shovel and place into a suitably labeled container for later disposal. **Do not dry sweep**

### 7. HANDLING AND STORAGE

Avoid generating dusts. Store out of direct sunlight. Store away from oxidizing agents. Keep containers closed, when not using the product. Store at 2-8°C and the reagent will be stable until the expiry date stated on the bottle and kit box labels. Store in original packages as approved by manufacturer.

### 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Exposure Standards**

No exposure standards are available for this product, however, the nature of the product is a dust therefore it is recommended that an exposure standard of 10 mg/m<sup>3</sup> for an eight hour exposure based upon Nuisance dust.

# Material Safety Data Sheet

Not classified as Hazardous according to criteria of Worksafe Australia  
Not classified as Hazardous according to Occupational Safety and Health Administration's (OSHA)  
Communication Standard, 29 CFR 1910.1200

## ALT (GPT)

### 8. EXPOSURE CONTROLS / PERSONAL PROTECTION (cont)

#### Engineering Controls

Maintain adequate ventilation at all times. In most circumstances natural ventilation systems are adequate .

#### Personal Protection Equipment

**GLOVES:** Not normally required, however, for people with sensitive skins the use of neoprene or nitrile is recommended.

**EYES:** Chemical goggles or spectacles to protect eyes.

**RESPIRATORY PROTECTION:** Avoid breathing of Dusts. Select and use respirators in accordance with AS/NZS 1715/1716. The use of a dust mask is not normally required, however, for high concentrations(for example in extremely dusty environments) use a half-face respirator fitted with a P1 filter. Filter capacity and respirator type depends on exposure levels and type of contaminant. If entering spaces where the airborne concentration of a contaminant is unknown then the use of a Self-contained breathing apparatus (SCBA) with positive pressure air supply complying with AS/NZS 1715 / 1716, or any other acceptable International Standard is recommended.

### 9. PHYSICAL AND CHEMICAL PROPERTIES

<b>Appearance:</b>	White to off-white powder with no odour.
<b>Boiling Point:</b>	Not available.
<b>Freezing Point:</b>	Not available.
<b>Vapour Pressure:</b>	Not available.
<b>Specific Gravity:</b>	Not available.
<b>Flash Point:</b>	Not applicable.
<b>Flammability Limits:</b>	Not applicable.
<b>Solubility in Water:</b>	Completely miscible.

#### Other Properties

**pH:** 7.65 - 7.95 at 19-22°C

### 10. STABILITY AND REACTIVITY

#### STABILITY:

Stable under normal conditions of use.

#### HAZARDOUS DECOMPOSITION PRODUCTS:

Decomposes on heating emitting oxides of carbon and noxious smoke.

#### HAZARDOUS POLYMERIZATION:

Will not occur.

#### INCOMPATIBILITIES:

Strong oxidizing agents.

#### CONDITIONS TO AVOID:

Generation of high dust levels and incompatibles.

### 11. TOXICOLOGICAL INFORMATION

There is no toxicological information available for this product, the following is expected:

Oral LD50(rat): > 5,000 mg/kg

Dermal LD50(rabbit): > 5,000 mg/kg

### 12. ECOLOGICAL INFORMATION

No ecological information available for this product. Do not dispose of large quantities to waterways, drains or sewers.

# Material Safety Data Sheet

Not classified as Hazardous according to criteria of Worksafe Australia  
Not classified as Hazardous according to Occupational Safety and Health Administration's (OSHA)  
Communication Standard, 29 CFR 1910.1200

## ALT (GPT)

### 13. DISPOSAL CONSIDERATIONS

Refer to appropriate authority in your State. Normally suitable for disposal by approved waste disposal agent.

### 14. TRANSPORT INFORMATION

**UN Number:** None allocated  
**Proper Shipping Name:** NONE ALLOCATED  
**Dangerous Goods Class:** None allocated  
**Subsidiary risk:** None allocated  
**Packing Group:** None allocated  
**Hazchem Code:** None allocated

Not classified as a Dangerous Good according to the Australian Code for the Transport of Dangerous Goods by Road and Rail (ADG Code) 6th Edition. Not classified as a Dangerous Good according to the UN, DOT(US), ICAO(IATA) or IMO(IMDG).

### 15. REGULATORY INFORMATION

**Poison Schedule:** None allocated

#### RISK PHRASES

None allocated

#### SAFETY PHRASES

S22 Do not breathe dust.

### 16. OTHER INFORMATION

#### Disclaimer

This MSDS summarises our best knowledge of the health and safety hazard information of the product and how to safely handle and use the product in the workplace. Each user should read this MSDS and consider the information in the context of how the product will be handled and used in the workplace including in conjunction with other products. If clarification or further information is needed to ensure that an appropriate risk assessment can be made, the user should contact this company.

Our responsibility for products sold is subject to our standard terms and conditions which are available on request.

© 2003 Thermo Electron Corporation. All rights reserved.

**Demand and Reply** are registered trademarks of Olympus Corporation, Lake Success, NY  
**All other trademarks are the property of Thermo Electron Corporation and its subsidiaries.**

License granted to make unlimited paper copies for internal use only.