

Material Safety Data Sheet

ThermoFisher
SCIENTIFIC

I CHEMICAL IDENTIFICATION		
NAME: Ammonia Standard Kit; Part Number: 1207		
DESCRIPTION: Ammonia solution		
Other Designations	Manufacturer	Emergency Contact
Ammonia Standards Ammonia Standard	NERL Diagnostics, LLC Part of Thermo Fisher Scientific Inc. 14 Almeida Avenue East Providence, RI 02914 Telephone: 1-800-556-7575 Fax: 1-401-438-2454	(401) 438-0386 (24 hour) CHEMTREC Phone: (800) 424-9300
II HAZARDOUS INGREDIENTS		III HEALTH HAZARDS
Acid		First Aid: Eye Contact: Flush with copious amounts of water. Get medical attention. Skin Contact: Remove contaminated clothing. Flush area with water. If irritation persists get medical attention. Ingestion: May be harmful if swallowed. If conscious, give large amounts of water. Seek medical attention. Inhalation: Do not breathe vapor.
IV SPECIAL PROTECTION INFORMATION		V SPECIAL PRECAUTIONS
Protective Clothing: Use appropriate eye/face protection and rubber gloves.		Wash exposed skin thoroughly with soap and water after handling. Do not reuse empty container, rinse container and put into trash.
VI SPILL OR LEAK PROCEDURES		VII REACTIVITY DATA
Absorb spilled material and dispose of in accordance with state and federal regulations. Waste disposal: Dispose of or incinerate in accordance with state and federal regulations.		Store at 2-8°C. Decomposes at temperatures above 8°C.
VIII FIRE AND EXPLOSION DATA		IX PHYSICAL DATA
Not a fire or explosion hazard.		
The information above is believed to be accurate and represents information currently available to us. NERL makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages from use of or reliance upon this information. This document is intended only as a guide to the appropriate precautionary handling of the material by properly trained persons.		