

# MATERIAL SAFETY DATA SHEET

Part of Thermo Fisher Scientific

## SECTION 1 - PRODUCT AND COMPANY IDENTIFICATION

Product Name:	Naturally occurring botanical materials		
Contents:	Naturally occurring botanical materials		
Company Identification	Microgenics Corporation	Company Phone Number	(510) 979-5000
Office Address	46360 Fremont Blvd.	Emergency Phone Number	(510) 979-5000
City, State, Zip Code	Fremont, CA 94538	CHEMTREC Phone Number, US	(800) 424-9300
Country	USA	CHEMTREC Phone Number, Europe	(202) 483-7616

## SECTION 2 - HAZARD IDENTIFICATION

CAUTION! May cause lung irritation if inhaled. By normal use the product is not regarded as hazardous to health. The product is not classified as flammable or harmful to the environment; may cause eye, skin, digestive track, and respiratory track irritation. The toxicological properties of this material have not been fully investigated. This product is not classified as dangerous in accordance with Directive 1999/45/EC or 67/548/EEC

## SECTION 3 - COMPOSITION/INFORMATION ON INGREDIENTS

Component: Name	CAS Number:	Percentage:	EINECS/ELINCS
Naturally occurring botanical materials	None	100	N/A

NOTE: Product contains one of the following: Bermuda Grass Smut Spores, Johnson Grass Smut Spores, Paper Mulberry Pollen, Ragweed Pollen, Bermuda Grass Pollen, Lycopodium Spores, Black Walnut Pollen, Pecan Pollen, Corn Pollen, and Walnut Shell Particles  
Component comments have no OSHA, Canadian, or EU hazardous ingredients above applicable threshold

## SECTION 4 - FIRST AID MEASURES

EYE CONTACT: Flush thoroughly with water for at least 15 minutes. If an irritation develops, seek medical advice  
SKIN CONTACT: In case of contact remove clothes, wash area with soap and flush with plenty of water. If an irritation develops, seek medical advice  
INHALATION: Remove to fresh air and seek medical advice  
INGESTION: Give one to two glasses of water and seek medical advice

## SECTION 5 - FIRE FIGHTING MEASURES

FLAMMABLE PROPERTIES: Flammable as dust. Finely dispersed particles form explosive mixtures in air.  
NFPA/HMIS Rating (estimated): Health: 1; Flammability: 0; Instability: 0  
PROTECTION FOR FIREFIGHTERS: Use extinguishing media appropriate for surrounding fire (water, foam, vaporizing liquid or multipurpose dry chemical)

## SECTION 6 - ACCIDENTAL RELEASE MEASURES

Areas covered with spilled microspheres may be slippery. If material is spilled or released, cordon off the area. Collect material by wiping the spill area with a paper towel or disposable wipe, and place materials into an appropriate container. When applicable use a vacuum cleaner equipped with a high efficiency particulate air (HEPA) filter and avoids inhaling fine particle dust. Discard collected material in containers suitable for proper disposal

## SECTION 7 - HANDLING AND STORAGE

Keep tightly sealed to prevent contamination; store at room temperature in a dry area. Avoid damaging or puncturing containers. When creating aerosols of fine particles use minimal and directional (away from the user) airflow.

## SECTION 8 - EXPOSURE CONTROL/PERSONAL PROTECTION

Occupational Exposure Limits "Particulates Not Otherwise Classified" (PNOC)

Chemical Name - Naturally occurring botanical materials

OSHA PEL	Cal/OSHA PEL	ACGIH TLV	AIHA	NIOSH REL
15mg/m <sup>3</sup> for total 5mg/m <sup>3</sup> for resp	10mg/m <sup>3</sup> - total dust 5 mg/m <sup>3</sup> - resp	10mg/m <sup>3</sup> for total 3mg/ m <sup>3</sup> for resp	N/A	N/A

EXPOSURE GUIDELINES: This product does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies. Safety glasses with side shields are suggested

ENGINEERING CONTROLS: Use minimal and directional (away from the user) airflow to minimize worker exposure

SKIN PROTECTION: Wear appropriate protective gloves to prevent skin exposure. No skin irritation expected

RESPIRATORY PROTECTION: Respiratory protection program that meets OSHA's 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149

GENERAL HYGIENE CONSIDERATIONS: Wear appropriate protective clothing to prevent skin exposure

## SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES

BOILING POINT:	N/A	SOLUBILITY IN WATER:	N/A	FLASH POINT:	N/A
MELTING POINT:	N/A	EVAPORATION RATE:	N/A	FLAMMABLE PROPERTIES:	N/A
PHYSICAL STATE:	solid	VAPOR PRESSURE:	N/A	AUTO IGNITION TEMP:	N/A
VAPOR DENSITY:	N/A	pH: N/A	Kow: N/A	DECOMPOSITION:	N/A
SPECIFIC GRAVITY:	N/A	APPEARANCE, COLOR, ODOR:		Yellow or brownish powder, faint floral odor	

**SECTION 10 - STABILITY AND REACTIVITY**

STABILITY: Stable at room temperature in closed containers under normal storage and handling conditions  
 INCOMPATIBILITY: None known  
 HAZARDOUS DECOMPOSITION PRODUCTS: None known  
 HAZARDOUS POLYMERIZATION: Will not occur

**SECTION 11 - TOXICOLOGICAL INFORMATION**

REPRODUCTIVE AND DEVELOPMENTAL TOXICITY: Material is not considered to be a reproductive or developmental toxicant  
 CARCINOGENICITY: Material is not listed by NTP, ACGIH, IARC or OSHA as a carcinogen. For more information refer to Section 2. Hazards identification

**SECTION 12 - ECOLOGICAL INFORMATION**

ECOTOXICITY Not regarded as dangerous for the environment LC/EC/IC50 > 100 mg/L  
 PERSISTENCE/DEGRADATION *No information found*  
 BIOACCUMULATION/ACCUMULATION *No information found*  
 MOBILITY IN ENVIRONMENTAL MEDIA Not soluble in water  
 ENVIRONMENTAL Not classified as harmful to the environment  
 OTHER ADVERSE EFFECTS *No information found*

**SECTION 13 - DISPOSAL CONSIDERATIONS**

Dispose of any waste residues according to prescribed federal, state, and local guidelines (e.g., to an appropriately permitted chemical waste incinerator)  
 None of the components are on RCRA P& RCRA U Series Wastes list. EWC waste code 18 01 05 discard as chemicals and medicines  
 Other Information: According to the European Waste Catalogue, Waste Codes are not product specific, but application specific. Waste codes should be assigned by the user based on the application for which the product was used

**SECTION 14 - TRANSPORTATION INFORMATION**

HAZARD CLASS: Not classified per US DOT, TDG, UNECE, ADR, ADN, IATA, ICAO, RID, IMDG/IMO, OTIF, Mexico

**SECTION 15 - REGULATORY INFORMATION**

TSCA 8(b)	SARA 302	SARA 313	CERCLA	RCRA 261.33	TSCA 8(d)	TSCA 12(b)	FIFRA	CWC	SARA 311/312	EPICRA	CAA- SNAP	FDA
no	no	no	no	no	no	no	no	no	no	no	no	no

Not subject to ITAR and EAR regulations

US State - N/A

California - Prop 65, Title 8 CCR, Title 22 CCR, ATHS, Occupational Exposure Limits, Ozone Depleting Compounds, Precursor Chemicals, SCAQMD, California No Significant Risk Level - *No information found*

None of the chemicals are on Prop 65 lists

Exposure Limits - *No information found*

European Union Directives - N/A

Canada Federal/Provincial Government - N/A

**WHMIS:**

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the MSDS contains all the information required by the CPR

**SECTION 16 -- OTHER INFORMATION**

Revision: 05

Date prepared: October 1, 2008

**ABBREVIATIONS**

EU: European Union  
 DOT: Department of Transportation  
 NIOSH: National Institute for Occupational Safety and Health  
 ACGIH: American Conference of Governmental Industrial Hygienists  
 IATA: International Air Transport Association  
 NTP: National Toxicology Program  
 OSHA: Occupational Safety and Health Administration  
 IARC: International Agency for Research on Cancer

*This MSDS is intended for review and guidance in the receipt, storage, handling, use and disposal of product purchased from Thermo Fisher Scientific, and for no other purpose. Use this product only as directed and in accordance with applicable instructions and warnings provided with the product. Please consult your institution's policies regarding use of this product. If you have obtained this MSDS other than in connection with the supply of this product from Thermo Fisher Scientific, this MSDS should be consulted for general information only, and should not be relied upon for any purpose. As with the use of all hazardous materials, you should in all instances follow the guidance of the MSDS provided or available with the specific product purchased.*

End of MSDS