

Material Safety Data Sheet

Enzymatic Creatinine Reagent

1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: Enzymatic Creatinine Reagent.

Catalog Numbers: TR35401

Use: This reagent comprises a two component system that is intended for in vitro quantitative determination of Creatinine in human serum, plasma or urine.

Thermo Fisher Scientific
189 - 199 Browns Rd
NOBLE PARK VIC 3174
AUSTRALIA

Tel: +61 3 9790 4100

Fax: +61 3 9790 4155

E-mail: info.clinicalchemistry@thermofisher.com

Contact Point

Australia

Quality Assurance Manager:

Tel: +61 3 9790 4100

Mon – Fri 9:00am to 5:00pm

Thermo Fisher Scientific

171 Industry Drive

Pittsburgh, PA, 15275

U.S.A

Tel: (800) 558 9115

Fax: (412) 788 6833

E-mail: info.clinicalchemistry@thermofisher.com

U.S.A

Chemtel

24 Hour Emergency Assistance

1-800-255-3924

813-248-0585

2. HAZARD IDENTIFICATION

Not classified as hazardous according to the ASCC/NOHSC/EU criteria.

Hazard Classification: NON-HAZARDOUS SUBSTANCE, NON DANGEROUS GOODS.

Hazard Category: None allocated.

RISK PHRASES

None allocated.

SAFETY PHRASES

S24/25 Avoid contact with skin and eyes.

Poison Schedule: None allocated [Aust].

3. COMPOSITION / INFORMATION ON INGREDIENTS

SUBSTANCE NAME

Proportion

CAS Number

All ingredients in this product are not hazardous according to the ASCC/NOHSC/EU criteria.

4. FIRST AID MEASURES

Swallowed: If swallowed, **DO NOT** induce vomiting. If conscious, give 1 to 2 glasses of water to drink. Seek immediate medical assistance.

Eye: If product enters the eyes, flush with plenty of water for at least 15 minutes, ensuring eye lids are held open. If irritation develops or persists, immediately transport to hospital or doctor.

Skin: If product contacts the skin, remove any contaminated clothing and wash skin thoroughly with soap and water. If irritation persists transport to hospital or doctor.

Inhaled: Move victim to fresh air. Apply resuscitation if victim is not breathing.

First Aid Facilities: Eye wash fountain, safety shower and normal wash room facilities.

Advice to Doctor: Treat symptomatically.

In case of poisoning, contact Poisons Information Centre

In Australia call Tel: 131126

In New Zealand Tel: 034747000

Material Safety Data Sheet

Enzymatic Creatinine Reagent

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media: Use extinguishing media suitable for surrounding fire situation.

Hazards from Combustion Products: Decomposes on heating emitting oxides of carbon, oxides of nitrogen, oxides of sulfur and noxious smoke.

Precautions for Fire Fighters and Special Protective Equipment: If safe to do so, move undamaged containers from fire area. Fire fighters to wear Self-contained breathing apparatus (SCBA) in confined spaces, in oxygen deficient atmospheres or if exposed to products of decomposition. Full protective clothing is also recommended.

Hazchem Code: None allocated.

Flammability: This material is not a combustible or flammable liquid.

6. ACCIDENTAL RELEASE MEASURES

Emergency Procedures

Keep unnecessary people away. Isolate hazard area and deny entry. If product spills onto floors it will represent a slip hazard, walk cautiously. Wear protective equipment to prevent skin and eye contact, as outlined under personal protection in this MSDS.

Methods and Materials for Containment and Cleanup Procedures

Dike area using with an absorbent such as diatomaceous earth - to prevent run off into drains and waterways. Throw further absorbent (diatomaceous earth or other inert material) on top of spill, then shovel up and seal in properly labeled containers for disposal.

7. HANDLING AND STORAGE

Precautions for Safe Handling

Do not swallow. Provide adequate ventilation. Avoid generating vapours.

Conditions for Safe Storage

Store in a cool place and out of direct sunlight. Store away from strong acids and oxidizers. **DO NOT FREEZE.** Light sensitive, store away from sunlight. Keep containers closed, when not using the product. When stored at 2 - 8°C the reagent will be stable until the expiry date on the bottle and kit box labels. Store in original packaging as approved by manufacturer.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure Standards

No exposure standards have been assigned by (NOHSC) for this product or any of the components.

Engineering Controls

Maintain adequate ventilation at all times. No other measures are required for this product.

Personal Protection Equipment

Gloves: Medical grade nitrile gloves are recommended for incidental handling.

Eyes: Chemical glasses or face shield for spills.

Respiratory Protection: Avoid breathing of vapours. The use of a respirator is not normally required, however, if entering spaces where the airborne concentration of a contaminant is unknown then the use of a Self-contained breathing apparatus (SCBA) with positive pressure air supply complying with AS/NZS 1715 / 1716, or any other acceptable International Standard is recommended. Select and use respirators in accordance with AS/NZS 1715/1716.

9. PHYSICAL AND CHEMICAL PROPERTIES

	Reagent 1	Reagent 2
Appearance:	Clear yellow coloured liquid with no odour	Clear yellow coloured liquid with no odour
Boiling Point:	Not available	Not available
Melting Point:	Not available	Not available
Vapour Pressure:	Not available	Not available
Specific Gravity:	Not available	Not available

Material Safety Data Sheet

Enzymatic Creatinine Reagent

9. PHYSICAL AND CHEMICAL PROPERTIES

Flash Point:	Not applicable	Not applicable
Flammability Limits:	Not applicable	Not applicable
Solubility in Water:	Completely miscible	Completely miscible

Other Properties	Reagent A	Reagent B
pH:	7.50 ± 0.1 @ 19 - 22°C	7.60 ± 0.1 @ 19 - 22°C

10. STABILITY AND REACTIVITY

Chemical Stability: Stable under normal conditions of use.

Conditions to Avoid: Eye and prolonged skin contact and mixing with incompatibles.

Incompatible Materials: Strong mineral acids and oxidizers.

Hazardous Decomposition Products: Decomposes on heating emitting oxides of carbon, oxides of nitrogen, oxides of sulfur and noxious smoke.

Hazardous Reactions: Will not occur.

11. TOXICOLOGICAL INFORMATION

There is no toxicological information available for this product, the following is expected:

Oral LD50(rat): > 5,000 mg/kg

Dermal LD50(rabbit): > 5,000 mg/kg

Acute Health Effects

Swallowed: Drinking large quantities of this product, may cause irritation to mouth, throat, which may lead to nausea, vomiting and diarrhoea.

Eye: May cause mild irritation to the eyes, with effects including: tearing and blurred vision. These effects are anticipated to be of a short acting nature and no long term injury is envisaged.

Skin: May cause mild irritation to the skin.

Inhaled: May cause mild irritation to the respiratory system.

Chronic: Chronic effects from repeated or long-term occupational exposure to this preparation are unknown. Prolonged or repeated skin contact may cause dermatitis.

12. ECOLOGICAL INFORMATION

No ecological information available for this product. Do not dispose of large quantities to waterways, drains or sewers.

13. DISPOSAL CONSIDERATIONS

Refer to appropriate authority in your State. Normally suitable for disposal by approved waste disposal agent.

14. TRANSPORT INFORMATION

UN Number: None allocated.

Proper Shipping Name: NONE ALLOCATED.

Dangerous Goods Class: None allocated.

Subsidiary risk: None allocated.

Packing Group: None allocated.

Hazchem Code: None allocated.

Road and Rail Transport:

Not classified as a Dangerous Good according to the United Nations Recommendations for the Transport of Dangerous Goods and Globally Harmonized System for the Classification and Labeling of Chemicals.

Material Safety Data Sheet

Enzymatic Creatinine Reagent

14. TRANSPORT INFORMATION

Air Transport:

Not classified as a Dangerous Good according to the International Civil Aviation Organization (ICAO) and International Air Transport Association (IATA) Dangerous Goods Regulations for transport by air.

Marine Transport:

Not classified as a Dangerous Good according to the International Maritime Organization Rules (Maritime Dangerous Goods Code - IMDG Code) for transport by sea.

15. REGULATORY INFORMATION

Poison Schedule: None allocated [Aust].

Inventory Status:

Australia (AICS)	Y
United States (TSCA)	Y
Canada (DSL)	Y
Europe (EINECS/ELINCS)	Y

Y = all ingredients are on the inventory.

If approved for European Communities use, this product is regulated under the In Vitro Diagnostic Medical Devices Directive (98/79/EC).

16. OTHER INFORMATION

Issue date: May, 2007.

Key Legend Information

NOHSC - National Occupational Health & Safety Commission [Aust]

SUSDP - Standard for the Uniform Scheduling of Drugs and Poisons {Poison Schedule} [Aust]

TWA - Time Weighted Average [Int]

STEL - Short Term Exposure Limit [Int]

AICS - Australian Inventory of Chemical Substances [Aust]

EPA - Environmental Protection Agency [Int]

NIOSH - National Institute for Occupational Safety and Health [US]

AS/NZS 1715 - Selection, Use and Maintenance of Respiratory Protective Devices. [Aust]

AS/NZS 1716 - Respiratory Protective Devices. [Aust]

Hazchem Code - Fire Fighter Designation [Aust]

IATA - International Aviation Transport Authority [Int]

ICAO - International Civil Aviation Organization [Int]

IMO - International Maritime Organisation. [Int]

IMDG - International Maritime Dangerous Goods [Int]

United Nations Recommendations for the Transport of Dangerous Goods and Globally Harmonized System for the Classification and Labeling of Chemicals. [Int]

EU - European Union

TSCA - Toxic Substances Control Act [US]

DSL - Domestic Chemical List [Can]

EINECS - European Inventory of Existing Commercial Chemical Substances [Int]

ELINCS - Existing List of Notified Chemical Substances. [Int]

[Aust] = Australia

[Int] = International

[US] = United States of America

[Can] = Canada

Principal References

Information supplied by manufacturer, reference sources including the public domain.

Material Safety Data Sheet

Enzymatic Creatinine Reagent

16. OTHER INFORMATION

Disclaimer

This MSDS summarises our best knowledge of the health and safety hazard information of the product and how to safely handle and use the product in the workplace. Each user should read this MSDS and consider the information in the context of how the product will be handled and used in the workplace including in conjunction with other products. If clarification or further information is needed to ensure that an appropriate risk assessment can be made, the user should contact this company.

Our responsibility for products sold is subject to our standard terms and conditions which are available on request.

© 2006-2007 Thermo Fisher Scientific. All rights reserved.

License granted to make unlimited paper copies for internal use only.

END OF MSDS