

Certificate of Analysis

Product Name: Aspire™ RP30 Desalting Tips

Product No.: 990-01

Lot No: 09050020

SPECIFICATIONS:

VISUAL: White to off-white SDVB resin housed in between two filters in a pipette tip

FRIT MATERIAL: Polyethylene (PE)

PIPETTE TIP MATERIAL: Polypropylene (PP)

BINDING CAPACITY (Saturated): Following the manufacturer's protocol to purify 60 µg of Bovine Serum Albumin (BSA) and recover at least 30 µg of BSA.

NUMBER OF TIPS/RACK: 96 tips per rack

RESULTS:

VISUAL: Pass

FRIT MATERIAL: Polyethylene (PE)

PIPETTE TIP MATERIAL: Polypropylene (PP)

BINDING CAPACITY: 41.89 µg Average

NUMBER OF TIPS/RACK: 96

Date Released: 01-26-2010

Expiration Date (16 month from release date): 07-26-2011

Quality Control Approved



Quality Manager