

Thermo Scientific HyClone Classical Media

Dulbecco's Modified Eagle's Medium/Nutrient Mixture F12 Ham's Liquid Media

Part Number	SH30023 SH30309		SH30261		SH30126	
Component Description	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
Inorganic Salts						
CALCIUM CHLORIDE (ANHY)	116.61	1.0507	116.61	1.0507	116.61	1.0507
COPPER SULFATE-5H2O	0.0013	0.00001	0.0013	0.00001	0.0013	0.00001
FERRIC NITRATE-9H2O	0.05	0.0001	0.05	0.0001	0.05	0.0001
FERROUS SULFATE-7H2O	0.417	0.0015	0.417	0.0015	0.417	0.0015
POTASSIUM CHLORIDE	311.8	4.1824	311.8	4.1824	311.8	4.1824
MAGNESIUM CHLORIDE (ANHY)	28.61	0.3005	0	0	28.61	0.3005
MAGNESIUM SULFATE (ANHY)	48.84	0.4056	84.95	0.7056	48.84	0.4056
SODIUM CHLORIDE	6999.5	119.7724	6999.5	119.7724	6999.5	119.7724
SODIUM PHOSPHATE DIBASIC, ANHY	71.02	0.5003	71.02	0.5003	71.02	0.5003
SODIUM PHOSPHATE MONOBASIC H2O	62.5	0.4529	62.5	0.4529	62.5	0.4529
ZINC SULFATE-7H2O	0.43	0.0015	0.43	0.0015	0.43	0.0015
Amino Acids						
L-ALANINE	4.45	0.05	4.45	0.05	4.45	0.05
L-ARGININE-HCL	147.5	0.7002	147.5	0.7002	147.5	0.7002
L-ASPARAGINE-H2O	7.501	0.05	7.501	0.05	7.501	0.05
L-ASPARTIC ACID	6.65	0.05	6.65	0.05	6.65	0.05
L-CYSTEINE-HCL-H2O	17.56	0.1	17.56	0.1	17.56	0.1
L-CYSTEINE-2HCL	31.29	0.0999	31.29	0.0999	31.29	0.0999
L-GLUTAMIC ACID	7.35	0.05	7.35	0.05	7.35	0.05
L-GLUTAMINE	365	2.4974	365	2.4974	0	0
GLYCINE	18.75	0.2498	18.75	0.2498	18.75	0.2498
L-HISTIDINE-HCL-H2O	31.48	0.1502	31.48	0.1502	31.48	0.1502
L-ISOLEUCINE	54.47	0.4153	54.47	0.4153	54.47	0.4153
L-LEUCINE	59.05	0.4502	59.05	0.4502	59.05	0.4502
L-LYSINE-HCL	91.26	0.4996	91.26	0.4996	91.26	0.4996
L-METHIONINE	17.24	0.1155	17.24	0.1155	17.24	0.1155
L-PHENYLALANINE	35.48	0.2148	35.48	0.2148	35.48	0.2148
L-PROLINE	17.25	0.1498	17.25	0.1498	17.25	0.1498
L-SERINE	26.25	0.2498	26.25	0.2498	26.25	0.2498
L-THREONINE	53.45	0.4487	53.45	0.4487	53.45	0.4487
L-TRYPTOPHAN	9.02	0.0442	9.02	0.0442	9.02	0.0442
L-TYROSINE-2NA-2H2O	55.79	0.2136	55.79	0.2136	55.79	0.2136
L-VALINE	52.85	0.4511	52.85	0.4511	52.85	0.4511
Vitamins						
D-BIOTIN	0.0036	0.00001	0.0036	0.00001	0.0036	0.00001
CALCIUM D-PANTOTHENATE	2.24	0.0047	2.24	0.0047	2.24	0.0047
CHOLINE CHLORIDE	8.98	0.0643	8.98	0.0643	8.98	0.0643
FOLIC ACID	2.65	0.006	2.65	0.006	2.65	0.006
MYO-INOSITOL	12.6	0.0699	12.6	0.0699	12.6	0.0699
NIACINAMIDE	2.02	0.0165	2.02	0.0165	2.02	0.0165
PYRIDOXAL HYDROCHLORIDE	0	0	0	0	2	0.0098
PYRIDOXINE-HCL	2.03	0.0099	2.03	0.0099	0.03	0.0002
RIBOFLAVIN	0.22	0.0006	0.22	0.0006	0.22	0.0006
THIAMINE-HCL	2.17	0.0064	2.17	0.0064	2.17	0.0064
VITAMIN B12	0.68	0.0005	0.68	0.0005	0.68	0.0005
D-GLUCOSE	3151	17.49	3151	17.49	3151	17.49
HEPES	3575	15.0021	3575	15.0021	3575	15.0021
HYPOXANTHINE-2NA	2.7	0.015	2.39	0.0132	2.7	0.015
LINOLEIC ACID	0.04	0.0001	0.04	0.0001	0.04	0.0001
DL-ALPHA-LIPOIC ACID	0.105	0.0005	0.105	0.0005	0.105	0.0005
PHENOL RED-NA	8.6	0.0229	8.1	0.0215	8.6	0.0229
PUTRESCINE-2HCL	0.08	0.0005	0.08	0.0005	0.08	0.0005
SODIUM PYRUVATE	55	0.4999	110	0.9996	55	0.4999
THYMIDINE	0.37	0.0015	0.37	0.0015	0.37	0.0015
Other						
SODIUM BICARBONATE	1200	14.284	2438	29.0204	1200	14.284

Thermo Scientific HyClone Classical Media

Dulbecco's Modified Eagle's Medium/Nutrient Mixture F12 Ham's Liquid Media

Part Number

SH30271

SH30272

	Component Description	mg/L	mmol/L	mg/L	mmol/L	
Inorganic Salts	CALCIUM CHLORIDE (ANHY)	116.61	1.0507	116.61	1.0507	
	COPPER SULFATE-5H2O	0.0013	0.00001	0.0013	0.00001	
	FERRIC NITRATE-9H2O	0.05	0.0001	0.05	0.0001	
	FERROUS SULFATE-7H2O	0.417	0.0015	0.417	0.0015	
	POTASSIUM CHLORIDE	311.8	4.1824	311.8	4.1824	
	MAGNESIUM CHLORIDE (ANHY)	28.61	0.3005	28.61	0.3005	
	MAGNESIUM SULFATE (ANHY)	48.84	0.4056	48.84	0.4056	
	SODIUM CHLORIDE	6999.5	119.7724	6999.5	119.7724	
	SODIUM PHOSPHATE DIBASIC, ANHY	71.02	0.5003	71.02	0.5003	
	SODIUM PHOSPHATE MONOBASIC H2O	62.5	0.4529	62.5	0.4529	
	ZINC SULFATE-7H2O	0.43	0.0015	0.43	0.0015	
	Amino Acids	L-ALANINE	4.45	0.05	4.45	0.05
		L-ARGININE-HCL	147.5	0.7002	147.5	0.7002
L-ASPARAGINE-H2O		7.501	0.05	7.501	0.05	
L-ASPARTIC ACID		6.65	0.05	6.65	0.05	
L-CYSTEINE-HCL-H2O		17.56	0.1	17.56	0.1	
L-CYSTINE-2HCL		31.29	0.0999	31.29	0.0999	
L-GLUTAMIC ACID		7.35	0.05	7.35	0.05	
L-GLUTAMINE		365	2.4974	365	2.4974	
GLYCINE		18.75	0.2498	18.75	0.2498	
L-HISTIDINE-HCL-H2O		31.48	0.1502	31.48	0.1502	
L-ISOLEUCINE		54.47	0.4153	54.47	0.4153	
L-LEUCINE		59.05	0.4502	59.05	0.4502	
L-LYSINE-HCL		91.26	0.4996	91.26	0.4996	
L-METHIONINE		17.24	0.1155	17.24	0.1155	
L-PHENYLALANINE		35.48	0.2148	35.48	0.2148	
L-PROLINE		17.25	0.1498	17.25	0.1498	
L-SERINE		26.25	0.2498	26.25	0.2498	
L-THREONINE		53.45	0.4487	53.45	0.4487	
L-TRYPTOPHAN		9.02	0.0442	9.02	0.0442	
L-TYROSINE-2NA-2H2O		55.79	0.2136	55.79	0.2136	
L-VALINE	52.85	0.4511	52.85	0.4511		
Vitamins	D-BIOTIN	0.0036	0.00001	0.0036	0.00001	
	CALCIUM D-PANTOTHENATE	2.24	0.0047	2.24	0.0047	
	CHOLINE CHLORIDE	8.98	0.0643	8.98	0.0643	
	FOLIC ACID	2.65	0.006	2.65	0.006	
	MYO-INOSITOL	12.6	0.0699	12.6	0.0699	
	NIACINAMIDE	2.02	0.0165	2.02	0.0165	
	PYRIDOXAL HYDROCHLORIDE	0	0	0	0	
	PYRIDOXINE-HCL	2.03	0.0099	2.03	0.0099	
	RIBOFLAVIN	0.22	0.0006	0.22	0.0006	
	THIAMINE-HCL	2.17	0.0064	2.17	0.0064	
	VITAMIN B12	0.68	0.0005	0.68	0.0005	
	D-GLUCOSE	3151	17.49	3151	17.49	
	HEPES	0	0	0	0	
	HYPOXANTHINE-2NA	2.7	0.015	2.7	0.015	
	LINOLEIC ACID	0.04	0.0001	0.04	0.0001	
	DL-ALPHA-LIPOIC ACID	0.105	0.0005	0.105	0.0005	
	PHENOL RED-NA	8.6	0.0229	0	0	
	PUTRESCINE-2HCL	0.08	0.0005	0.08	0.0005	
	SODIUM PYRUVATE	55	0.4999	55	0.4999	
	THYMIDINE	0.37	0.0015	0.37	0.0015	
	Other	SODIUM BICARBONATE	1200	14.284	1200	14.284

Cell Culture & BioProcessing

925 West 1800 South
Logan, UT
84321

In Americas/Asia
435-792-8000
435-792-8001 fax

In Europe
+32 53 85 71 80
+32 53 85 74 31 fax

www.thermo.com/hyclone

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific