

Thermo Scientific HyClone Classical Media

Phosphate Buffered Saline Solutions

Part Number	SH30013		SH30028		SH30264	
Component Description	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
CALCIUM CHLORIDE (ANHY)	0	0	0	0	100	0.9011
POTASSIUM CHLORIDE	200	2.6828	200	2.6828	200	2.6828
POTASSIUM PHOSPHATE MONOBASIC ANHY	200	1.4696	200	1.4696	200	1.4696
MAGNESIUM CHLORIDE-6H2O	0	0	0	0	100	0.4919
SODIUM CHLORIDE	8000	136.8925	8000	136.8925	8000	136.8925
SODIUM PHOSPHATE DIBASIC, ANHY	1150	8.1009	1150	8.1009	1150	8.1009
SODIUM PHOSPHATE DIBASIC 7H2O	0	0	0	0	0	0

Phosphate Buffered Saline Solutions

Part Number	SH30378		SH30597		SH30256	
Component Description	mg/L	mmol/L	mg/L	mmol/L	mg/L	mmol/L
CALCIUM CHLORIDE (ANHY)	0	0	1000	9.0106	0	0
POTASSIUM CHLORIDE	2000	26.8276	2000	26.8276	0	0
POTASSIUM PHOSPHATE MONOBASIC ANHY	2000	14.6962	2000	14.6962	144	1.0581
MAGNESIUM CHLORIDE-6H2O	0	0	1000	4.9188	0	0
SODIUM CHLORIDE	80000	1368.9254	80000	1368.9254	9000	154.0041
SODIUM PHOSPHATE DIBASIC, ANHY	0	0	11500	81.0087	795	5.6002
SODIUM PHOSPHATE DIBASIC 7H2O	21600	80.576	0	0	0	0

Phosphate Buffered Saline Solutions

Part Number	SH30258	
Component Description	mg/L	mmol/L
CALCIUM CHLORIDE (ANHY)	0	0
POTASSIUM CHLORIDE	0	0
POTASSIUM PHOSPHATE MONOBASIC ANHY	1440	10.5812
MAGNESIUM CHLORIDE-6H2O	0	0
SODIUM CHLORIDE	90000	1540.0411
SODIUM PHOSPHATE DIBASIC, ANHY	7950	56.0017
SODIUM PHOSPHATE DIBASIC 7H2O	0	0