

	<p><b>SV-288</b>                  Aluminum Rotor                  8-Place, 36 ml                  Tube Dimensions: 25 × 90 mm                  Vertical: 0°                  Maximum Speed: 20,000 rpm                  K Factor: 210                  Rotor Mass: 8.3 kg</p>		RCF (g)	R (cm)
		Max.	40,301	9.02
		Min.	28,908	6.47

		RC-28S	RC-26 Plus	RC-24 <sup>1</sup>	RC-5C Plus RC-5C <sup>1</sup>	RC-5B Plus RC-5B <sup>1</sup> RC-5 <sup>1</sup>	RC-2B <sup>1</sup> RC-2 <sup>1</sup>	SS-3 <sup>1</sup> SS-4 <sup>1</sup>	RC-3B <sup>1</sup> RC-3B Plus RC-3C <sup>1</sup> RC-3C Plus
<b>SV-288</b>	Max Speed (rpm)	20,000	20,000	20,000	20,000	20,000	n/a	n/a	n/a
	Max RCF (g)	40,301	40,301	40,301	40,301	40,301	n/a	n/a	n/a

Basic Catalog No. 09203	Complete Catalog No. 09201	Accessories
Includes: SV-288 Rotor Lubricant <b>61556</b> Rotor Extractor <b>01031</b> Instruction Manual <b>09057</b>	Includes Basic Catalog <b>09203</b> Plus: ULTRACRIMP® Tubes, 33 ml, Polyallomer (8 pkgs, 25/pkg) <b>03899</b> <sup>2</sup> Rotor Cap (8 pkgs, 1/pkg) <b>03995</b> Tube Racks (1 pkg, 2/pkg) <b>03763</b> ULTRACRIMP® Tube Removal Tool <b>12826</b>	ULTRACRIMP® Sealing Tool and Crimp Gauge <b>03920</b> Crimp Gauge, Replacement <b>03919</b> ULTRACRIMP® Plugs and Caps, Extra (25 each/pkg) <b>03999</b> Tool Kit for ULTRACRIMP® <b>12770</b> . Includes: Rotor Holding Fixture <b>52540</b> , Torque Wrench <b>52580</b> , Socket <b>52882</b> Cap Wrench (for Polyallomer Thin-Walled Tubes) <b>09054</b> Tube Trimmer, optional <b>09073</b> Rotor Holding Fixture <b>52540</b> <sup>3</sup>

**B-2**

**TUBES**

CAPACITY (ml)		DESCRIPTION	MAX SPEED (rpm)	DIMENSIONS (mm)	TUBE CAT. NO. (Qty/Pkg)	SEALING ASSEMBLY*		ADAPTER REQ'D		NOTES
NOM-INAL	FILL					CAT. NO. (Qty/Pkg)	DESCRIPTION	CAT. NO. (Qty/Pkg)	NO. OF PLACES	
36	33	Polyallomer Thin-Walled Tube	20,000	25 × 87	<b>03141</b> (25)	<b>09082</b> (1) <b>52482</b> (8)	Body Cap Plug	— —	— —	Requires Body Cap <b>09082</b> , Tube Plug <b>52482</b> , Rotor Holding Fixture <b>52540</b> , and Wrench <b>09054</b> . Tube Trimmer <b>09073</b> is optional.
33	33	Polyallomer ULTRACRIMP® Tube	20,000	25 × 90	<b>03899</b> (25)	Included with Tubes	Plug and Aluminum Cap	—	—	Requires Rotor Cap <b>03995</b> , Tube Racks <b>03763</b> , ULTRACRIMP® Sealing Tool <b>03920</b> and Tool Kit <b>12770</b> .
20	17	Polyallomer Thin-Walled Tube	20,000	25 × 51	<b>03288</b> (25)	<b>09082</b> (1) <b>52482</b> (8)	Body Cap Plug	<b>00421</b> (1)	1	Requires Body Cap <b>09082</b> , Tube Plug <b>52482</b> , Rotor Holding Fixture <b>52540</b> , and Wrench <b>09054</b> .

<sup>(1)</sup>Instrument no longer available.

<sup>(2)</sup>Tubes supplied with equal number of plugs and caps.

<sup>(3)</sup>Included in Tool Kit **12770** for ULTRACRIMP® tubes. Must be ordered separately for use with tubes **03288** and **03141**.

\*Sealing Assembly replacement parts may be available. See section E for availability.